Peel and Place Dressing Kit

Safety Information and Application Instructions

Only for use with 3M™ V.A.C.® Therapy Systems

Rx Only

Table of Contents

3M™ V.A.C.® Therapy System Safety Information5

 Indications for Use5

 Contraindications6

 Warnings.....6

 Additional Warnings for the 3M™ V.A.C.® Peel and Place Dressing9

 Precautions10

 Considerations for Transitioning 3M™ V.A.C.® Therapy Into Home Care.....12

3M™ V.A.C.® Peel and Place Dressing Specific Application Instructions.....13

 Component Identification.....13

 3M™ V.A.C.® Peel and Place Dressing Change Interval14

 3M™ V.A.C.® Peel and Place Dressing Removal.....14

 3M™ V.A.C.® Peel and Place Dressing Specific Information.....15

 3M™ V.A.C.® Peel and Place Dressing Application Instructions16

 3M™ SensaT.R.A.C.™ Pad Application.....19

 3M™ V.A.C.® Therapy Application with the 3M™ V.A.C.® Peel and Place Dressing.....20

 Symbols Used.....22

3M™ V.A.C.® Therapy System Safety Information

Disposable components of the 3M™ V.A.C.® Therapy System are to be used as indicated on the associated product labeling. 3M™ V.A.C.® Therapy Unit canisters are packaged sterile or fluid path sterile. All disposable components of the 3M™ V.A.C.® Therapy System are for single use only. To help ensure safe and effective use, the 3M™ V.A.C.® Peel and Place Dressing is to be used only with 3M™ V.A.C.® Therapy Units.

Re-use of disposable components may result in wound contamination, infection and / or failure of the wound to heal.

The decision to use clean versus sterile / aseptic technique is dependent upon wound pathophysiology, physician / clinician preference and institutional protocol.

Important: As with any prescription medical device, failure to consult a physician and carefully read and follow all therapy unit and dressing instructions and safety information prior to each use may lead to improper product performance and the potential for serious or fatal injury. Do not adjust therapy unit settings or perform therapy application without directions from / or supervision by the treating physician.

Indications for Use

The 3M™ V.A.C.® Peel and Place Dressing is intended to be applied in an acute, extended, or home care setting where product application is conducted by or under the supervision of a qualified healthcare professional. It is an accessory to the 3M™ Negative Pressure Wound Therapy System.

When used on open wounds, the system is intended to create an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudate and infectious material. Open wound types include: chronic, acute, traumatic, subacute and dehiscent wounds, partial-thickness burns, ulcers (such as diabetic, pressure or venous insufficiency), flaps and grafts.

Contraindications

- Do not place foam dressings of the 3M™ V.A.C.® Therapy System directly in contact with exposed blood vessels, anastomotic sites, organs or nerves.

Note: Refer to Warnings section for additional information concerning Bleeding.

- 3M™ V.A.C.® Therapy is contraindicated for patients with:

- Malignancy in the wound
- Untreated osteomyelitis

Note: Refer to Warnings section for Osteomyelitis information.

- Non-enteric and unexplored fistulas
- Necrotic tissue with eschar present

Note: After debridement of necrotic tissue and complete removal of eschar, 3M™ V.A.C.® Therapy may be used.

Warnings

Bleeding: With or without using 3M™ V.A.C.® Therapy, certain patients are at high risk of bleeding complications. The following types of patients are at increased risk of bleeding, which, if uncontrolled, could be potentially fatal:

- Patients who have weakened or friable blood vessels or organs in or around the wound as a result of, but not limited to:
 - Suturing of the blood vessel (native anastomoses or grafts) / organ
 - Infection
 - Trauma
 - Radiation
- Patients without adequate wound hemostasis
- Patients who have been administered anticoagulants or platelet aggregation inhibitors
- Patients who do not have adequate tissue coverage over vascular structures

If 3M™ V.A.C.® Therapy is prescribed for patients who have an increased risk of bleeding complications, they should be treated and monitored in a care setting deemed appropriate by the treating physician.

If active bleeding develops suddenly or in large amounts during 3M™ V.A.C.® Therapy, or if frank (bright red) blood is seen in the tubing or in the canister, immediately stop 3M™ V.A.C.® Therapy, leave dressing in place, take measures to stop the bleeding and seek immediate medical assistance. The 3M™ V.A.C.® Therapy Units and dressings should not be used to prevent, minimize or stop vascular bleeding.

- **Protect Vessels and Organs:** All exposed or superficial vessels and organs in or around the wound must be completely covered and protected prior to the administration of 3M™ V.A.C.® Therapy. Always ensure that 3M™ V.A.C.® Foam Dressings do not come in direct contact with vessels or organs. Use of a thick layer of natural tissue should provide the most effective protection. If a thick layer of natural tissue is not available or is not surgically possible, multiple layers of meshed non-adherent material or bio-engineered tissue may be considered as an alternative, if deemed by the treating physician to provide a complete protective barrier. If using non-adherent materials, ensure they are secured in a manner that will maintain their protective position throughout therapy. Consideration should also be given to the negative pressure wound therapy setting and therapy mode used when initiating therapy.

Caution should be taken when treating large wounds that may contain hidden vessels which may not be readily apparent. The patient should be closely monitored for bleeding in a care setting deemed appropriate by the treating physician.

- **Infected Blood Vessels:** Infection may erode blood vessels and weaken the vascular wall which may increase susceptibility to vessel damage through abrasion or manipulation. **Infected blood vessels are at risk of complications, including bleeding, which, if uncontrolled, could be potentially fatal. Extreme caution should be used when 3M™ V.A.C.® Therapy is applied in close proximity to infected or potentially infected blood vessels.** (Refer to Protect Vessels and Organs section above). The patient should be closely monitored for bleeding in a care setting deemed appropriate by the treating physician.
- **Hemostasis, Anticoagulants, and Platelet Aggregation Inhibitors:** Patients without adequate wound hemostasis have an increased risk of bleeding, which, if uncontrolled, could be potentially fatal. These patients should be treated and monitored in a care setting deemed appropriate by the treating physician.
Caution should be used in treating patients on doses of anticoagulants or platelet aggregation inhibitors thought to increase their risk for bleeding (relative to the type and complexity of the wound). Consideration should be given to the negative pressure wound therapy setting and therapy mode used when initiating therapy.
- **Hemostatic Agents Applied at the Wound Site:** Non-sutured hemostatic agents (for example, bone wax, absorbable gelatin sponge or spray wound sealant) may, if disrupted, increase the risk of bleeding, which, if uncontrolled, could be potentially fatal. Protect against dislodging such agents. Consideration should be given to the negative pressure wound therapy setting and therapy mode used when initiating therapy.
- **Sharp Edges:** Bone fragments or sharp edges could puncture protective barriers, vessels or organs, causing injury. Any injury could cause bleeding, which, if uncontrolled, could be potentially fatal. Beware of possible shifting in the relative position of tissues, vessels or organs within the wound that might increase the possibility of contact with sharp edges. Sharp edges or bone fragments must be covered or eliminated from the wound area, to prevent them from puncturing blood vessels or organs before the application of 3M™ V.A.C.® Therapy. Where possible, completely smooth and cover any residual edges to decrease the risk of serious or fatal injury, should shifting of structures occur. Use caution when removing dressing components from the wound so that wound tissue is not damaged by unprotected sharp edges.

1000 mL Canister and the Risk of Excessive Fluid Loss, Including Blood: Consider the size and weight of the patient, patient condition (patients with a high risk of bleeding or on patients unable to tolerate loss of fluid volume, including children or the elderly), wound type, monitoring capability and care setting when using the 1000 mL canister. Patients should be closely monitored for excessive fluid loss and dehydration, as well as frank blood in the canister. The 1000 mL canister is recommended for acute care use only.

- **3M™ V.A.C.® Rx4 Therapy Unit:** The 3M™ V.A.C.® Rx4 Therapy Unit provides four independent therapy channels that may accommodate either 500 or 1000 mL canisters. **When using multiple channels on multiple wounds, DO NOT USE the 1000 mL canister for patients at high risk for excessive fluid loss.**

Infected Wounds: Infected wounds should be monitored closely and may require more frequent dressing changes than non-infected wounds, dependent upon factors such as wound conditions and treatment goals. Refer to application instructions for details regarding dressing change frequency. As with any wound treatment, clinicians and patients / caregivers should frequently monitor the patient's wound, periwound tissue and exudate for signs of infection, worsening infection or other complications. Some signs of infection are fever, tenderness, redness, swelling, itching, rash, increased warmth in the wound or periwound area, purulent discharge or strong odor. Infection can be serious, and can lead to complications such as pain, discomfort, fever, gangrene, toxic shock, septic shock and / or fatal injury. Some signs or complications of systemic infection are nausea, vomiting, diarrhea, headache, dizziness, fainting, sore throat with swelling of the mucus membranes, disorientation, high fever, refractory and / or orthostatic hypotension or erythroderma (a sunburn-like rash). **If there are any signs of the onset of systemic infection or advancing infection at the wound site, contact the treating physician immediately to determine if 3M™ V.A.C.® Therapy should be discontinued.** For wound infections relating to blood vessels, please also refer to the section titled **Infected Blood Vessels**.

Osteomyelitis: 3M™ V.A.C.® Therapy should NOT be initiated on a wound with untreated osteomyelitis. Consideration should be given to thorough debridement of all necrotic, non-viable tissue, including infected bone (if necessary) and appropriate antibiotic therapy.

Dressing Placement: Always use 3M™ V.A.C.® Dressings from sterile packages that have not been opened or damaged. Do not place any portion of the dressing into blind / unexplored tunnels. Do not force dressings into any area of the wound, as this may damage tissue, alter the delivery of negative pressure wound therapy or hinder exudate and foam removal.

3M™ V.A.C.® Foam Dressings are radiolucent, not detectable on X-Ray.

Dressing Removal: Gently remove the dressing from the wound. Regardless of treatment modality, disruption of the new granulation tissue during any dressing change may result in bleeding at the wound site. Minor bleeding may be observed and considered expected. However, patients with increased risk of bleeding, as described in the **Bleeding** section, have a potential for more serious bleeding from the wound site. **If significant bleeding develops, immediately discontinue the use of the 3M™ V.A.C.® Therapy System, take measures to stop the bleeding, and do not remove the foam dressing until the treating clinician is consulted. Do not resume the use of the 3M™ V.A.C.® Therapy System until adequate hemostasis has been achieved and the patient is not at risk for continued bleeding.**

Keep 3M™ V.A.C.® Therapy On: Never leave a 3M™ V.A.C.® Dressing in place without active 3M™ V.A.C.® Therapy for more than two hours. If therapy is off for more than two hours, remove the old dressing and irrigate the wound. Either apply a new 3M™ V.A.C.® Dressing from an unopened sterile package and restart 3M™ V.A.C.® Therapy, or apply an alternative dressing at the direction of the treating physician.

Acrylic Adhesive and Silicone Layer: The 3M™ V.A.C.® Peel and Place Dressing has an acrylic adhesive coating and a silicone layer, which may present a risk of an adverse reaction in patients who are allergic or hypersensitive to acrylic adhesives or silicone. If a patient has a known allergy or hypersensitivity to such materials, do not use the 3M™ V.A.C.® Therapy System. If any signs of allergic reaction or hypersensitivity develop, such as redness, swelling, rash, urticaria or significant pruritus, discontinue use and consult a physician immediately. If bronchospasm or more serious signs of allergic reaction appear, seek immediate medical assistance.

Defibrillation: Remove the 3M™ V.A.C.® Dressing if defibrillation is required in the area of dressing placement. Failure to remove the dressing may inhibit transmission of electrical energy and / or patient resuscitation.

Magnetic Resonance Imaging (MRI) - 3M™ V.A.C.® Therapy Unit: The 3M™ V.A.C.® Therapy Unit is MR unsafe. Do not take the 3M™ V.A.C.® Therapy Unit into the MR environment.

Magnetic Resonance Imaging (MRI) - 3M™ V.A.C.® Dressings: 3M™ V.A.C.® Dressings can typically remain on the patient with minimal risk in an MR environment, assuming that use of the 3M™ V.A.C.® Therapy System is not interrupted for more than two hours (refer to Keep 3M™ V.A.C.® Therapy On section).

Hyperbaric Oxygen Therapy (HBO): Do not take the 3M™ V.A.C.® Therapy Unit into a hyperbaric oxygen chamber. The 3M™ V.A.C.® Therapy Unit is not designed for this environment and **should be considered a fire hazard**. After disconnecting the 3M™ V.A.C.® Therapy Unit, either (i) replace the 3M™ V.A.C.® Dressing with another HBO compatible material during the hyperbaric treatment or (ii) cover the unclamped end of the 3M™ V.A.C.® Tubing with dry gauze. For HBO therapy, the 3M™ V.A.C.® Tubing must not be clamped. Never leave a 3M™ V.A.C.® Dressing in place without active 3M™ V.A.C.® Therapy for more than two hours (refer to Keep 3M™ V.A.C.® Therapy On section).

Note: *The 3M™ V.A.C.® Granufoam™ Bridge Dressing contains additional synthetic materials which may pose a risk during HBO Therapy.*

Additional Warnings for the 3M™ V.A.C.® Peel and Place Dressing

Wound Fillers: DO NOT use any wound fillers with the 3M™ V.A.C.® Peel and Place Dressing. If a wound filler is required, select an alternative dressing to use.

Cutting the Foam: Do not cut the foam section of the dressing.

Wound Depth: Do not use the 3M™ V.A.C.® Peel and Place Small Dressing on wounds with a greater depth than 2 cm. Do not use the 3M™ V.A.C.® Peel and Place Medium Dressing on wounds with a greater depth than 4 cm. Do not use the 3M™ V.A.C.® Peel and Place Large Dressing on wounds with a greater depth than 6 cm. If the wound is deeper than 6 cm, use an alternate 3M™ Dressing.

Tunnels and Undermining: Do not place the 3M™ V.A.C.® Peel and Place Dressing into a wound with any tunnels or with undermining greater than 2 cm. If clinician considers wound to be inappropriate for the 3M™ V.A.C.® Peel and Place Dressing, please use an alternate 3M™ Dressing.

Pressure Settings and Continuous Therapy Only: The 3M™ V.A.C.® Peel and Place Dressing must only be used at pressure settings between 75 and 150 mmHg in Continuous mode.

Instillation / 3M™ Veraflo™ Therapy: DO NOT use with 3M™ Veraflo™ Therapy (Instillation) provided by the 3M™ V.A.C.® Ulta Therapy Unit. Instillation into the wound dressed with the 3M™ V.A.C.® Peel and Place Dressing may result in fluid leaks which may result in maceration.

Open Abdomen: DO NOT USE the 3M™ V.A.C.® Peel and Place Dressing over the open abdomen or with 3M™ AbThera™ Open Abdomen Negative Pressure Therapy. Use over the open abdomen may result in an inability to maintain a negative pressure seal.

Use over Umbilicus: If the 3M™ V.A.C.® Peel and Place Dressing covers the umbilicus, the umbilicus should be filled with an antimicrobial petroleum gauze prior to dressing application.

Precautions

Standard Precautions: To reduce the risk of transmission of bloodborne pathogens, apply standard precautions for infection control with all patients, per institutional protocol, regardless of their diagnosis or presumed infection status. In addition to gloves, use gown and goggles if exposure to body fluids is likely.

Continuous versus Intermittent / 3M™ Dynamic Pressure Control™ Therapy: Continuous rather than intermittent / 3M™ Dynamic Pressure Control™ Therapy is required.

Patient Size and Weight: The size and weight of the patient should be considered when prescribing 3M™ V.A.C.® Therapy. Infants, children, certain small adults and elderly patients should be closely monitored for fluid loss and dehydration. Also, patients with highly exuding wounds or large wounds in relation to the patient size and weight should be closely monitored, as these patients have a risk of excessive fluid loss and dehydration. When monitoring fluid output, consider the volume of fluid in both the tubing and canister.

Spinal Cord Injury: In the event a patient experiences autonomic dysreflexia (sudden changes in blood pressure or heart rate in response to stimulation of the sympathetic nervous system), discontinue 3M™ V.A.C.® Therapy to help minimize sensory stimulation and seek immediate medical assistance.

Bradycardia: To minimize the risk of bradycardia, 3M™ V.A.C.® Therapy must not be placed in proximity to the vagus nerve.

Enteric Fistulas: Wounds with enteric fistulas require special precautions to optimize 3M™ V.A.C.® Therapy. 3M™ V.A.C.® Therapy is not recommended if enteric fistula effluent management or containment is the sole goal of therapy.

Protect Periwound Skin: The 3M™ V.A.C.® Peel and Place Dressing incorporates a hybrid acrylic and silicone drape so use of a skin preparation product is not necessary.

- Multiple layers of drape may decrease the moisture vapor transmission rate, which may increase the risk of maceration.
- If any signs of irritation or sensitivity to the drape, foam, or tubing assembly appear, discontinue use and consult a physician.
- To avoid trauma to the periwound skin, do not pull or stretch the drape over the foam dressing during drape application.

Circumferential Dressing Application: Avoid use of circumferential dressings. Extreme care should be taken not to stretch or pull the drape when securing it, but let it attach loosely and stabilize the edges with an elastic wrap, if necessary. When using circumferential drape applications, it is crucial to systematically and recurrently palpate distal pulses and assess distal circulatory status. If circulatory compromise is suspected, discontinue therapy, remove dressing, and contact a treating physician.

3M™ V.A.C.® Therapy Unit Pressure Excursions: In rare instances, tubing blockages with the 3M™ V.A.C.® Therapy Unit may result in brief vacuum excursions to more than 250 mmHg negative pressure. Resolve alarm conditions immediately. Refer to the Therapy Unit User Guide or Manual or contact your 3M representative for additional information.

In addition to these general warnings and precautions for 3M™ V.A.C.® Therapy, additional warnings and precautions apply to certain 3M™ V.A.C.® specialty dressings and 3M™ V.A.C.® Therapy Units. Please refer to the specific product instructions for use and labeling prior to application for complete safety information, dressing application instructions, specific therapy settings, and the procedure for connection to the 3M™ V.A.C.® Therapy Unit.

Additional Precautions for the 3M™ V.A.C.® Rx4 Therapy Unit

- **Single Patient Use:** 3M™ V.A.C.® Rx4 Therapy Unit is not intended for use on multiple patients simultaneously as this may pose additional risks of cross contamination at the dressing site.
- **1000 mL Canisters and the Risk of Excessive Fluid Loss, Including Blood:** When using multiple channels on multiple wounds, DO NOT USE the 1000 mL canister for patients at high risk for excessive fluid loss.
- **Not for Home Use:** The 3M™ V.A.C.® Rx4 Therapy Unit is intended for use in acute care settings and other professional healthcare environments where product use is conducted by or under the supervision of a qualified healthcare professional.

Considerations for Transitioning 3M™ V.A.C.® Therapy Into Home Care

Warning: Patients with an increased risk of bleeding complications should be treated and monitored in a care setting deemed appropriate by the treating physician.

- **The 3M™ V.A.C.® Rx4 and 3M™ V.A.C.® Ulta Therapy Systems are NOT intended for home use.**
- If there is a need to continue 3M™ V.A.C.® Therapy when a patient transitions home, consider using one of the 3M Therapy Systems approved for the post-acute environment, such as:
 - 3M™ ActiV.A.C.™ Therapy Unit
 - 3M™ V.A.C.® Simplicity Therapy Unit

Refer to the safety information included with those devices for important information.


In addition to the contraindications, warnings and precautions for use of 3M™ V.A.C.® Therapy, consider the following before prescribing 3M™ V.A.C.® Therapy for use in the home care setting.

- **The Patient's Situation:**
 - Clinical condition (adequate hemostasis and a low risk of active and / or large amounts of bleeding at the wound site)
 - Home environment (patient or family member / caregiver able to read and understand safety labeling, able to respond to alarms, able to follow instructions for use)
- **The Patient's Wound:**
 - Must be assessed for exposed vessels, anastomotic sites, organs, and nerves. Adequate protection must be present.
- **The 3M™ V.A.C.® Therapy System Canister Size:**
 - The 1000 mL canister is NOT intended for use in the home.
- **Labeling:**
 - The prescribing physician and health care clinician should be familiar with the 3M™ V.A.C.® Therapy instructional materials that accompany the therapy unit into the home. The prescribing physician and / or healthcare clinician should carefully review these materials with the patient and patient's caregiver.
 - 3M offers in-service and training programs for use of 3M™ V.A.C.® Therapy. Contact your local 3M representative. In the U.S., call 1-800-275-4524 for scheduling.

For additional and most current information, please see 3M's website at HCBGRegulatory.3M.com.

3M™ V.A.C.® Peel and Place Dressing Specific Application Instructions

Component Identification

3M™ V.A.C.® Peel and Place Dressing

3M™ SensaT.R.A.C.™ Pad

Disposable components of the 3M™ V.A.C.® Therapy System including the 3M™ V.A.C.® Peel and Place Dressing and tubing are packaged sterile. 3M™ V.A.C.® Therapy Unit canisters are packaged sterile or fluid path sterile. All disposable components of the 3M™ V.A.C.® Therapy System are for single use only. To help ensure safe and effective use, the 3M™ V.A.C.® Peel and Place Dressing is to be used only with 3M™ V.A.C.® Therapy Units.

The decision to use clean versus sterile/aseptic technique is dependent upon wound pathophysiology, clinician preference and institutional protocol.

Always consult a physician and review and follow 3M™ V.A.C.® Therapy Safety Information and 3M™ V.A.C.® Therapy Unit instructions.

Devices for use in conjunction with the 3M™ V.A.C.® Peel and Place Dressing

For **acute care settings**, the 3M™ V.A.C.® Peel and Place Dressing may be used with:

- 3M™ V.A.C.® Ulta Therapy Unit (3M™ V.A.C.® Therapy only)
- 3M™ V.A.C.® Rx4 Therapy Unit
- 3M™ ActiV.A.C.™ Therapy Unit
- 3M™ V.A.C.® Simplicity Therapy Unit

For **home care settings**, the 3M™ V.A.C.® Peel and Place Dressing may be used only with:

- 3M™ ActiV.A.C.™ Therapy Unit
- 3M™ V.A.C.® Simplicity Therapy Unit

3M™ V.A.C.® Peel and Place Dressing Change Interval

Wounds being treated with the 3M™ V.A.C.® Therapy System should be monitored on a regular basis. In a monitored, non-infected wound, **3M™ V.A.C.® Peel and Place Dressings may be left in place for up to 7 days with the frequency adjusted by the clinician as appropriate.** Infected wounds must be monitored often and very closely. For these wounds, dressings may need to be changed more often; the dressing change intervals should be based on a continuing evaluation of wound condition and the patient's clinical presentation, rather than a fixed schedule.

Warning: Review all 3M™ V.A.C.® Therapy System Safety Information before beginning Wound Preparation.

1. Remove and discard previous dressing per institution protocol. Thoroughly inspect wound to ensure all pieces of dressing components have been removed.

3M™ V.A.C.® Peel and Place Dressing Removal

2. Gently remove an existing 3M™ V.A.C.® Peel and Place Dressing according to the following procedure:
 - a. Raise the tubing connectors above the level of the therapy unit.
 - b. Close clamp on the dressing tubing.
 - c. Disconnect canister tubing from dressing tubing.
 - d. Allow the therapy unit to pull the exudate in the canister tube into the canister, then close the clamp on the canister tubing.
 - e. Press THERAPY ON/OFF to deactivate the 3M™ V.A.C.® Therapy Unit. Wait 15-30 seconds to allow foam to decompress.
 - f. Gently remove the 3M™ V.A.C.® Peel and Place Dressing from the wound.
 - g. Discard disposables according to institutional or state regulations.

Note: If dressing adheres to wound, consider introducing sterile water or normal saline into the dressing, waiting 15-30 minutes, then gently remove the dressing from the wound.

3. Debride all necrotic, non-viable tissue, including bone, eschar or hardened slough, as prescribed by physician.
4. Perform thorough wound and periwound area cleaning per physician order or institution protocol prior to each dressing application.
5. Ensure adequate hemostasis has been achieved (refer to Warnings, Bleeding section, Hemostasis, Anticoagulants and Platelet Aggregation Inhibitors).
6. Prior to dressing placement, protect vessels and organs (refer to Warnings, Bleeding section, Protect Vessels and Organs).
7. Sharp edges or bone fragments must be eliminated from wound area or covered (refer to Warnings, Bleeding section, Sharp Edges).



3M™ V.A.C.® Peel and Place Dressing Specific Information

- The 3M™ V.A.C.® Peel and Place Dressing has a protective interface layer so the foam area can be placed over both the wound and surrounding intact skin. The interface layer allows the dressing to be used for up to 7 days.
- The drape border can be lifted and repositioned once during initial application to simplify dressing placement.
- Refer to the **3M™ V.A.C.® Peel and Place Dressing Application Instructions section** for appropriate **wound types** and **depths**. **An alternative 3M dressing may be required for some wounds.**
- Refer to the **3M™ V.A.C.® Therapy Application** with 3M™ V.A.C.® Peel and Place Dressing section for **recommended pressure settings**. Negative pressure must be sufficient to stretch the thin foam into the wound bed.
- The 3M™ V.A.C.® Peel and Place Dressing is intended to be used **without sizing the drape border** for easy application, however, **trimming may be appropriate for some applications**, at the discretion of the clinician. **Leave adequate drape border for sealing.**

Warning: Do not cut the foam section of the dressing.

3M™ V.A.C.® Peel and Place Dressing Application Instructions

1. **Assess wound** for any **undermining** or **tunneling**.


Warning: Do not place the 3M™ V.A.C.® Peel and Place Dressing into a wound with any tunnels or with undermining greater than 2 cm. If clinician considers wound inappropriate for the 3M™ V.A.C.® Peel and Place Dressing, please **use an alternate 3M dressing**.





2. Assess wound dimensions and depth.

- **First**, choose dressing size based on **complete coverage** of the wound bed.


Note: The 3M™ V.A.C.® Peel and Place Dressing foam should **extend beyond the edge of the wound** to ensure full coverage after 3M™ V.A.C.® Therapy is applied.

- **Second**, check **maximum wound depth** recommended for the chosen 3M™ V.A.C.® Peel and Place Dressing size. **If the wound is deeper than 6 cm, use an alternate 3M dressing.**


3. Remove the release liners to expose adhesive.


4. Place the dressing onto the wound with the pre-cut hole facing up and positioned for the desired location of the 3M™ SensaT.R.A.C.™ Pad. The hole does not need to be at the center of the wound.

Note: Do not use wound fillers.

5. **Check for complete coverage of the wound.** The foam edges will be pulled inward when negative pressure is applied.

Note: *If the dressing exposes wound edges after initial application, lift and readjust the dressing. **Do not lift and reapply during therapy.***

6. **Assess dressing placement.** If needed, the dressing **may be lifted**, repositioned and reapplied once for a better fit.

Note: *Superficial or retention sutures should be covered with a single layer of non-adherent material placed between the sutures and the drape.*

7. **Smooth** any wrinkles or creases in the dressing to prevent leaks.


Note: *Drape can be peeled back and reapplied during initial placement to address affected part of the dressing.*

8. While holding down the edge of the dressing, remove the perforated handling bars and pat down to ensure an occlusive seal.

3M™ SensaT.R.A.C.™ Pad Application


Note: Do not cut off the pad or insert the tubing directly into the foam dressing. This may occlude the tubing and cause the 3M™ V.A.C.® Therapy Unit to alarm.


1. The 3M™ V.A.C.® Peel and Place Dressing comes pre-cut with a hole for the 3M™ SensaT.R.A.C.™ Pad.

Note: Give particular consideration to tubing positioning to allow for optimal flow and avoid placement over bony prominences or within creases in the tissue.

2. Apply pad, which has a central disc and a surrounding outer adhesive skirt.


- a. Remove both backing Layers 1 and 2 to expose adhesive.
- b. Place pad opening in central disc directly over hole in drape.
- c. Apply gentle pressure on the central disc and outer skirt to ensure complete adhesion of the pad.
- d. Pull back on blue tab to remove pad stabilization layer.

3M™ V.A.C.® Therapy Application with the 3M™ V.A.C.® Peel and Place Dressing

Warning: Review all 3M™ V.A.C.® Therapy System Safety Information before initiating 3M™ V.A.C.® Therapy.

1. Remove 3M™ V.A.C.® Canister from packaging and insert into the 3M™ V.A.C.® Therapy Unit until it locks into place.

Note: If the canister is not fully engaged, the 3M™ V.A.C.® Therapy Unit will alarm.

2. Connect 3M™ SensaT.R.A.C.™ Pad tubing to canister tubing and ensure clamp on each tube is open. Position clamps away from patient.
3. Turn on power to the 3M™ V.A.C.® Therapy Unit and select the prescribed therapy setting.

Warning: The 3M™ V.A.C.® Peel and Place Dressing **must only be used at pressure settings between 75 and 150 mmHg in Continuous mode.**

4. Initiate 3M™ V.A.C.® Therapy. Assess dressing to ensure seal integrity. The dressing should be collapsed and have a wrinkled appearance. There should be no hissing sounds. For the 3M™ ActiV.A.C.™, 3M™ V.A.C.® Ultra and 3M™ V.A.C.® Rx4 Therapy Systems, use the 3M™ Seal Check™ Feature to verify that the rate of air leakage is below the alarm threshold. If there is any evidence of non-integrity, check 3M™ SensaT.R.A.C.™ Pad and drape seals, tubing connections and canister insertion, and ensure clamps are open. If there is a leak, pat down the dressing or lift and reapply to ensure that there are no creases in the drape.

Note: Do not lift and reapply drape edges during therapy (only during initial placement).

5. Secure excess tubing to prevent interference with patient mobility.

Note: Refer to unit specific user guide or manual and / or quick reference guide for information regarding alarms.

Note: If a leak source is identified, patch with additional drape to ensure seal integrity.


Note: If the wound is over a bony prominence or in areas where weight bearing may exert additional pressure or stress to the underlying tissues, a pressure redistribution (pressure relief) surface or device should be used to optimize patient offloading.

Warning: Never leave a 3M™ V.A.C.® Dressing in place without active 3M™ V.A.C.® Therapy for more than two hours. If therapy is off for more than two hours, remove the old dressing and irrigate the wound. Either apply a new 3M™ V.A.C.® Dressing from an unopened sterile package and restart 3M™ V.A.C.® Therapy, or apply an alternative dressing such as wet to moist gauze, as approved during times of extreme need by treating clinician.

Symbols Used

Symbol Title	Symbol	Description and Reference
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Container quantity		Indicates quantity of medical devices in container
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Do not re-sterilize		Indicates a medical device that is not to be re-sterilized. Source: ISO 15223, 5.2.6
Do not use if package is damaged and consult instructions for use		Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. Source: ISO 15223, 5.2.8
Keep dry		Indicates a medical device that needs to be protected from moisture. Source: ISO 15223, 5.3.4
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Rx Only		Indicates the U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Refer to instruction manual / booklet		To signify that the instruction manual / booklet must be read. ISO 7010-M002

* This symbol is blue on packaging.

Single sterile barrier system		Indicates a single sterile barrier system. Source: ISO 15223, 5.2.11
Sterilized using ethylene oxide		Indicates a medical device that has been sterilized using ethylene oxide. Source: ISO 15223, 5.2.3
Sterilized using irradiation		Indicates a medical device that has been sterilized using irradiation. Source: ISO 15223, 5.2.4
Unique device identifier		Indicates a carrier that contains Unique Device Identifier information. Source: ISO 15223, 5.7.10
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
KCI now part of 3M		

For more information see, [HCBGregulatory.3M.com](https://www.hcbgregulatory.3m.com)

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**KCI USA, Inc.**