PREVENA PLUS[™] INCISION MANAGEMENT SYSTEM

PREVENA PLUS[™] CUSTOMIZABLE[™] DRESSING WITH SENSAT.R.A.C.[™] TECHNOLOGY

INSTRUCTIONS FOR USE FOR CLINICIANS ONLY

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

KCI HEREBY DISCLAIMS ALL EXPRESS OR IMPLIED WARRANTIES, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE KCI PRODUCT(S) DESCRIBED IN THIS PUBLICATION. ANY WRITTEN WARRANTY OFFERED BY KCI SHALL BE EXPRESSLY SET FORTH IN THIS PUBLICATION OR INCLUDED WITH THE PRODUCT. UNDER NO CIRCUMSTANCES SHALL KCI BE LIABLE FOR ANY INDIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES AND EXPENSES, INCLUDING DAMAGES OR INJURY TO PERSON OR PROPERTY, DUE IN WHOLE OR IN PART TO THE USE OF THE PRODUCT OTHER THAN THOSE FOR WHICH DISCLAIMER OF WARRANTY OR LIMITATION OF LIABILITY IS EXPRESSLY PROHIBITED BY SPECIFIC, APPLICABLE LAW. NO PERSON HAS THE AUTHORITY TO BIND KCI TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH IN THIS PARAGRAPH.

Descriptions or specifications in KCI printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties except as set forth in the written limited warranty included with this product. Information in this publication may be subject to change at any time. Contact KCI for updates.

TABLE OF CONTENTS

Product Description and Indication for Use	5
Important Information for Users.....	5
Optimum Use Conditions.....	6
Contraindication.....	6
Warnings.....	6
Precautions.....	7
PREVENA PLUS™ Incision Management System Site Preparation	8
Drain Tubes and Pain Management Control Devices.....	8
PREVENA PLUS™ Incision Management System with PREVENA PLUS™ CUSTOMIZABLE™ Dressing.....	9
PREVENA PLUS™ CUSTOMIZABLE™ Dressing Applications.....	10
Dressing Application for Linear Incisions	11
Dressing Application for Non-linear Incisions.....	13
Dressing Application for Intersecting Incisions.....	14
Drape Application	15
SENSAT.R.A.C.™ Pad Application.....	17
PREVENA PLUS™ Canister Installation.....	18
Connecting the Dressing to the PREVENA PLUS™ 125 Therapy Unit.....	18
Beginning Therapy.....	19
Unit Troubleshooting	20
Therapy Life Indicators.....	20
Duration of PREVENA PLUS™ Therapy	21
Alerts.....	21
Correcting a Leak Condition.....	23
Check Canister Tubing Connection.....	23
Indications That a Leak Condition Has Been Corrected.....	24
Battery Charging.....	24
Canister Removal and Replacement	26
Using the PREVENA PLUS™ Dressing with KCI V.A.C.® Therapy Units.....	26
Connecting the PREVENA PLUS™ Dressing to V.A.C.® Therapy Units.....	26
Setting Negative Pressure on the V.A.C.® Therapy Units	27
Alarm Resolutions	27
PREVENA PLUS™ 125 Therapy Unit Disposal	27
Instructions for Patients	27
Daily Use	27
Sleeping	27
Showering and Bathing	28

Strenuous Activity	28
Cleaning.....	28
Dressing Removal.....	28
Warnings and Important Information for Users - PREVENA PLUS™ 125 Therapy Unit	29
Customer Contact Information	30
PREVENA PLUS™ Therapy System Specifications	30
PREVENA PLUS™ 125 Therapy Unit Electromagnetic Compatibility	31
Included Power Supplies	36
Bibliography of Published Studies	37
Symbols Used	495

INSTRUCTIONS FOR USE PREVENA PLUS™ INCISION MANAGEMENT SYSTEM

KCI CUSTOMER CONTACT INFORMATION IS LOCATED IN THE BACK OF THIS GUIDE.

PRODUCT DESCRIPTION AND INDICATION FOR USE

The PREVENA PLUS™ Incision Management System is intended to manage the environment of closed surgical incisions and surrounding intact skin in patients at risk for developing post-operative complications, such as infection, by maintaining a closed environment via the application of a negative pressure wound therapy system to the incision. The PREVENA™ Dressing skin interface layer with silver reduces microbial colonization in the fabric.

The system consists of:

A PREVENA PLUS™ Dressing and a source of negative pressure, which may be one of the following KCI therapy units:

- PREVENA PLUS™ 125 Therapy Unit
- ACTIV.A.C.™ Therapy Unit
- V.A.C.ULTA™ Therapy Unit

The PREVENA™ Dressing is only intended to be used with negative pressure and should not be used as a stand-alone device.

Clinical studies have been conducted on KCI Negative Pressure Incision Management Systems. Refer to the **Bibliography of Published Studies** in the back of this guide.

IMPORTANT INFORMATION FOR USERS

WARNING: DO NOT use with V.A.C. VERAFLU™ Therapy (Instillation) provided by the V.A.C.ULTA™ Therapy Unit. Instillation into the incision site may result in pooling of fluid which may result in maceration.

WARNING: The V.A.C.ULTA™ Therapy Unit is only indicated for use in acute care settings. Before transitioning the patient to home care, this therapy unit must be replaced with one for home use, such as the PREVENA PLUS™ 125 or ACTIV.A.C.™ Therapy Unit.

For pressure settings and connection information for use of the PREVENA PLUS™ Dressings with the V.A.C.® Therapy Units listed above, see the **Using the PREVENA PLUS™ Dressing with KCI V.A.C.® Therapy Units** section.

CAUTION: The PREVENA PLUS™ Incision Management System should be applied and removed only by qualified physicians or nurses.

As with any prescription medical device, failure to carefully read and follow all instructions and safety information prior to use may lead to improper product performance.

The PREVENA PLUS™ Incision Management System dressings and therapy unit canisters are disposable and are for single use only. Re-use of disposable components may result in wound contamination, infection and/or failure of the wound to heal.

OPTIMUM USE CONDITIONS

For maximum benefit the PREVENA PLUS™ Incision Management System should be applied immediately post surgery to clean surgically closed wounds. It is to be continuously applied for a minimum of two days up to a maximum of seven days. It can transition home with the patient.

The PREVENA PLUS™ Incision Management System will not be effective in addressing complications associated with:

- ischemia to the incision or incision area
- untreated or inadequately treated infection
- inadequate hemostasis of the incision
- cellulitis of the incision area

The PREVENA PLUS™ Incision Management System should not be used to treat open or dehiscenced surgical wounds.

The V.A.C.® Therapy System should be considered for treatment of these wounds.

Consider using the smallest available canister for the selected V.A.C.® Therapy Unit.

The PREVENA PLUS™ Incision Management System should be used with caution in the following patients:

- patients with fragile skin surrounding the incision as they may experience skin or tissue damage upon removal of the PREVENA PLUS™ Dressing
- patients who are at an increased risk of bleeding from the incision associated with the use of anticoagulants and/or platelet aggregation inhibitors

CONTRAINDICATION

- sensitivity to silver

WARNINGS

The PREVENA PLUS™ Incision Management System is not intended to manage open or dehiscenced wounds.

DO NOT use with V.A.C. VERAFLU™ Therapy (Instillation) provided by the V.A.C. ULTA™ Therapy Unit. Instillation into the incision site may result in pooling of fluid which may result in maceration.

Bleeding: Before applying the PREVENA PLUS™ Incision Management System to patients who are at risk of bleeding complications due to the operative procedure or concomitant therapies and/or co-morbidities, ensure that hemostasis has been achieved and all tissue planes have been approximated. If active bleeding develops suddenly or in large amounts during therapy, or if frank blood is seen in the tubing or in the canister, the patient should leave the PREVENA PLUS™ Dressing in place, turn off the therapy unit and seek immediate emergency medical assistance.

Infected Wounds: As with any wound treatment, clinicians and patients/caregivers should frequently monitor the patient's wound, periwound tissue and exudate for signs of infection or other complications. Some signs of infection are fever, tenderness, redness, swelling, itching, rash, increased warmth in the wound or periwound area, purulent discharge or strong odor. Infection can be serious, and can lead to complications such as pain, discomfort, fever, gangrene, toxic shock, septic shock and/or fatal injury. Some signs or complications of systemic infection are nausea, vomiting, diarrhea,

headache, dizziness, fainting, sore throat with swelling of the mucus membranes, disorientation, high fever, refractory and/or orthostatic hypotension or erythroderma (a sunburn-like rash). Silver in the interface layer of the PREVENA PLUS™ Dressing is not intended to treat infection, but to reduce bacterial colonization in the fabric. **If infection develops, PREVENA PLUS™ Therapy should be discontinued until the infection is treated.**

Allergic Response: The PREVENA PLUS™ Dressing has an acrylic adhesive coating and a skin interface layer with silver, which may present a risk of an adverse reaction in patients who are allergic or hypersensitive to acrylic adhesives or silver. If a patient has a known allergy or hypersensitivity to these materials, do not use the PREVENA PLUS™ Dressings. If any signs of allergic reaction, irritation or hypersensitivity develop, such as redness, swelling, rash, urticaria or significant pruritus, patient should consult a physician immediately. If bronchospasm or more serious signs of allergic reaction appear, the patient should turn off the therapy unit and seek immediate emergency medical assistance.

Defibrillation: Remove the PREVENA PLUS™ Dressing if defibrillation is required in the area of dressing placement. Failure to remove the dressing may inhibit transmission of electrical energy and/or patient resuscitation.

Magnetic Resonance Imaging (MRI): All KCI Therapy Units, including the PREVENA PLUS™ 125 Therapy Unit, are MR unsafe. Do not take therapy units into the MR environment. The PREVENA PLUS™ Dressings can typically remain on the patient with minimal risk in an MR environment. Interruption of PREVENA PLUS™ Therapy during MRI may reduce the effectiveness of the PREVENA PLUS™ Incision Management System. The PREVENA PLUS™ Dressings pose no known hazards in an MR environment with the following conditions of use: static magnetic field of 3 Tesla or less, spatial gradient field of 720 Gauss/cm or less and maximum whole-body-averaged specific absorption rate (SAR) of 3W/kg for 15 minutes of scanning.

Diagnostic Imaging: The PREVENA PLUS™ Dressing contains metallic silver that may impair visualization with certain imaging modalities.

Hyperbaric Oxygen Therapy (HBO): Do not take the PREVENA PLUS™ 125 Therapy Unit, V.A.C.® Therapy Units or PREVENA PLUS™ Dressings into a hyperbaric oxygen chamber. They are not designed for this environment and **should be considered a fire hazard**. If PREVENA PLUS™ Therapy is reinitiated after HBO treatment, do not readhere the same dressing; a new dressing must be applied.

Canister Full: If at any time while using the PREVENA PLUS™ Incision Management System the canister becomes full of fluid, indicated by a therapy unit alert or visual inspection, the patient should turn off the therapy unit and contact the treating physician for additional instruction.

Standard Operation: Do not use accessories or materials not provided with the PREVENA PLUS™ Incision Management System. For a list of acceptable therapy units with which PREVENA PLUS™ Dressings may be used, see the **Product Description and Indications for Use** section.

PRECAUTIONS

Standard Precautions: To reduce the risk of transmission of bloodborne pathogens, apply standard precautions for infection control with all patients, per institutional protocol, regardless of their diagnosis or presumed infection status.

Circumferential Dressing Application: Avoid applying the PREVENA PLUS™ Dressing circumferentially. In cases where the clinician determines that the benefits of applying the PREVENA PLUS™ Dressing circumferentially outweigh the risk of circulatory compromise, extreme care should be taken not to stretch or pull the dressing when securing it. Attach the dressing loosely and stabilize

edges with an elastic wrap if necessary. It is crucial to systematically and recurrently palpate distal pulses and assess distal circulatory status. If circulatory compromise is suspected, discontinue therapy and remove dressing.

Electrodes or Conductive Gel: Do not allow the PREVENA PLUS™ Dressing to come in contact with EKG or other electrodes or conductive gels during electronic monitoring or when taking electronic measurements.

Dressing Components: The PREVENA PLUS™ Dressing contains ionic silver (0.019%). Application of products containing silver may cause temporary tissue discoloration.

- Always use PREVENA PLUS™ Dressings and canisters from sterile packages that have not been opened or damaged.
- All dressing components and canisters of the PREVENA PLUS™ Incision Management System are for single use only. Do not re-use any component of this system.
- To avoid trauma to the skin, do not pull or stretch the adhesive border of the dressing during application.

Compressive Garments or Dressings: Avoid tight compressive garments or dressings (such as surgical bras, elastic bandage wraps or abdominal binders) to prevent forcibly pressing the PREVENA PLUS™ Dressing into soft tissue.

PREVENA PLUS™ INCISION MANAGEMENT SYSTEM SITE PREPARATION

1. Prior to surgery, shave or clip the surgical area where the dressing will be applied to improve dressing adhesion and seal integrity.
2. Gather all items needed for application:
 - sterile wound cleaning solution, e.g. water, saline or alcohol
 - sterile gauze or other material to clean application site
 - all components of the PREVENA PLUS™ Incision Management System (dressing and therapy unit)
3. After surgery, cleanse the application site with sterile gauze and sterile wound cleaning solution using a circular motion beginning at the center of the surgical area and extending outward to ensure that the site is free of foreign material.
4. Pat the application site dry with sterile gauze. To ensure proper adhesion, the application site must be completely dry before dressing is applied.


DRAIN TUBES AND PAIN MANAGEMENT CONTROL DEVICES

The PREVENA PLUS™ Incision Management System can be used with both drain tubes and pain devices, provided the dressing is not placed over tubing where it exits the skin. Surgical drains must be routed under the skin beyond the boundary of the dressing and function independently of the PREVENA PLUS™ Incision Management System.

NOTE: While the concomitant use of surgical drains is allowable with the PREVENA PLUS™ Incision Management System, the system must not be used as an outlet or reservoir for the drain.

PREVENA PLUS™ INCISION MANAGEMENT SYSTEM WITH PREVENA PLUS™ CUSTOMIZABLE™ DRESSING


The PREVENA PLUS™ CUSTOMIZABLE™ System can be used for linear/non-linear incision applications and contains the following single use, disposable components.

PREVENA PLUS™ CUSTOMIZABLE™ Dressing (folded) - a specially designed dressing for application to the surgical area


Sealing Strips - used to create a continuous adhesive seal around the PREVENA PLUS™ CUSTOMIZABLE™ Dressing


KCI Drapes - used to cover dressing and provide a seal for negative pressure

Ruler - used to measure incision before cutting. The removeable label may be used as needed to record date of dressing application or removal.

SENSAT.R.A.C.™ Pad - used to connect the PREVENA PLUS™ CUSTOMIZABLE™ Dressing to a KCI therapy unit


PREVENA PLUS™ 125 Therapy Unit - delivers negative pressure to the surgical area. The unit is battery and electrically powered. The non-sterile PREVENA PLUS™ 125 Therapy Unit Carrying Case is provided to facilitate patient mobility.

PREVENA PLUS™ 150 ml Canister - a sterile reservoir for collection of wound fluids

PREVENA PLUS™ Connector - Connects the SENSAT.R.A.C.™ Pad tubing to the canister


PREVENA PLUS™ 125 Therapy Unit Power Supply with Power Cord - a charging system provided with the PREVENA PLUS™ 125 Therapy Unit to charge the internal battery.

PREVENA PLUS™ CUSTOMIZABLE™ DRESSING APPLICATIONS


The PREVENA PLUS™ CUSTOMIZABLE™ Dressing can be used for the following wound types:

- linear incisions
- non-linear incisions
- intersecting incisions


The following figures show completed dressing applications using the PREVENA PLUS™ CUSTOMIZABLE™ Dressing, KCI Drape and SENSAT.R.A.C.™ Pad.

Linear Incision
Dressing Application

Non-linear Incision
Dressing Application

Intersecting Incision
Dressing Application



WARNING: DO NOT use with V.A.C. VERAFLU™ Therapy (Instillation) provided by the V.A.C. ULTA™ Therapy Unit. Instillation into the incision site may result in pooling of fluid which may result in maceration.

DRESSING APPLICATION FOR LINEAR INCISIONS


CAUTION: If the dressing covers the umbilicus, the umbilicus must first be fully filled with an anti-microbial petroleum gauze prior to dressing application.

1. Open the sterile dressing package, remove and unfold foam dressing using aseptic technique. Do not use if package has been torn or the sterile seal has been compromised.
2. Use the enclosed ruler to measure incision before cutting.
3. Cut foam dressing to the length necessary to cover the incision.

CAUTION: Do not cut the foam over the wound, as fragments may fall into the wound. Away from wound site, rub foam edges to remove any fragments or loose particles that may fall into or be left in the wound upon dressing removal.



4. Remove bottom release liner 1 on the foam dressing to expose adhesive seal and white fabric located on the bottom of the foam dressing. The white fabric covers the incision.


5. Add sealing strips to the ends of the foam dressing.


NOTE: Use the PREVENA PLUS™ CUSTOMIZABLE™ Dressing to size the sealing strip as needed to fit the width of the dressing, approximately 3 in (8 cm).

- Liner 1 **must be removed from the dressing** before applying sealing strips.
- Remove one of the liners from the sealing strip. Sealing strips are tacky on both sides.
- Place sealing strip on end of foam dressing (creating a continuous seal). The strip should be placed half on dressing and half off dressing and overlap the adhesive seal on each side of the dressing. This ensures a seal all around the edge of the dressing.
- Remove the remaining liner from the sealing strip. Repeat strip application on other end of foam dressing.



6. With sealing strips applied, position the foam dressing until the **adhesive seal and white fabric are facing down**.


7. Center the foam dressing over the incision/surgical closure. Ensure the **adhesive seal and white fabric are facing down**.

8. Press down around the length of the adhesive seal to ensure adhesion to the skin.

9. Peel off release liner 2 from the top of the adhesive seals.





10. Once the foam dressing has been applied to the incision, continue with the drape application (see **Drape Application** section).

DRESSING APPLICATION FOR NON-LINEAR INCISIONS

CAUTION: If the dressing covers the umbilicus, the umbilicus must first be fully filled with an anti-microbial petroleum gauze prior to dressing application.

1. Open the sterile dressing package, remove and unfold foam dressing using aseptic technique. Do not use if package has been torn or the sterile seal has been compromised.
2. Use enclosed ruler to measure incision before cutting.
3. Cut foam dressing to the length and shape necessary to cover the incision.

CAUTION: Do not cut the foam over the wound, as fragments may fall into the wound. Away from wound site, rub foam edges to remove any fragments or loose particles that may fall into or be left in the wound upon dressing removal.

4. Make straight cuts at the desired angles to ensure the foam dressing pieces and seals make a flush/even connection with **no gaps**.

CAUTION: Gaps at foam junctions may potentially cause erythema due to pinching of the skin when negative pressure is applied.


5. Use a sealing strip to bridge where the foam dressing pieces meet. Sealing strip should have equal overlap on both foam dressing pieces.

NOTE: Use the PREVENA PLUS™ CUSTOMIZABLE™ Dressing to size the sealing strip as needed to fit the width of the dressing, approximately 3 in (8 cm).

- Liner 1 **must be removed from the dressing** before applying sealing strips.
- Place sealing strip on end of foam dressing (creating a continuous seal). The strip should be placed half on dressing and half off dressing and overlap the adhesive seal on each side of the dressing. This ensures a seal all around the edge of the dressing.

NOTE: Having exposed adhesive material around all edges of the dressing provides a greater seal when the drape is applied.

- Remove the remaining liner from the sealing strip. Repeat strip application on other end of foam dressing.

6. With sealing strips applied, position the foam dressing until the **adhesive seal and white fabric are facing down**.


NOTE: It is acceptable if sealing strips cover the incision line.

7. Place the foam dressing over the incision/surgical closure. Ensure the **adhesive seal and white fabric are facing down**.
8. After foam dressing placement is completed, ensure all liners are removed before continuing with the drape application (see **Drape Application** section).

DRESSING APPLICATION FOR INTERSECTING INCISIONS

CAUTION: If the dressing covers the umbilicus, the umbilicus must first be fully filled with an anti-microbial petroleum gauze prior to dressing application.

1. Open the sterile dressing package, remove and unfold foam dressing using aseptic technique. Do not use if package has been torn or the sterile seal has been compromised.
2. Use enclosed ruler to measure incision before cutting.



3. Cut dressing as needed to fit incision.

CAUTION: Do not cut the foam over the wound, as fragments may fall into the wound. Away from wound site, rub foam edges to remove any fragments or loose particles that may fall into or be left in the wound upon dressing removal.

4. Ensure the foam dressings and seals meet with **no gaps**.

CAUTION: Gaps at foam junctions may potentially cause erythema due to pinching of the skin when negative pressure is applied.

CAUTION: Do not overlap foam dressing pieces on top of each other.

5. Use sealing strips at the ends of the foam dressing to create a continuous seal. Sealing strip should have equal overlap on both foam dressing pieces.

NOTE: Use the PREVENA PLUS™ CUSTOMIZABLE™ Dressing to size the sealing strip as needed to fit the width of the dressing, approximately 3 in (8 cm).

- Liner 1 **must be removed from the dressing** before applying sealing strips.

6. With sealing strips applied, position the foam dressing until the **adhesive seal and white fabric are facing down**.


NOTE: It is acceptable if sealing strips cover the incision line.

7. Place the foam dressing over the incision/surgical closure. Ensure the **adhesive seal and white fabric are facing down**.
8. After foam dressing placement is completed, ensure all liners are removed before continuing with the drape application (see **Drape Application** section).


DRAPE APPLICATION

CAUTION: Do not cover the umbilicus with drape. Before or after placement, ensure a small hole is cut in the drape in relation to the umbilicus.


Multiple pieces of drape may be needed when using PREVENA PLUS™ CUSTOMIZABLE™ Dressings depending on incision size. Drapes should always overlap each other by at least 1.2 in (3 cm). All dressing liners must be removed before any drape is used.


1. Trim KCI Drape as necessary. Ensure a minimum overlap of 1.2 in (3 cm) when using multiple pieces of KCI Drape. There is no drape under the blue handling bars.

2. Carefully remove the middle section of layer 1 from KCI Drape to expose adhesive. Leave the two end sections of layer 1 intact to maintain drape stability during application. The KCI Drape should be held by the blue handling bars.


3. With adhesive side down, apply KCI Drape over foam and onto intact skin. Ensure KCI Drape covers a minimum of 1.2 in (3 cm) of intact periwound skin.


4. Lift the drape at both ends and remove the two remaining end sections of layer 1. Smooth out any wrinkles in the KCI Drape.


CAUTION: Do not press blue handling bars down onto the sealing material.

5. Carefully remove striped layer 2 from KCI Drape.

6. Remove the two blue handling bars.


7. Add additional drape as needed to cover foam dressings in various configurations, as shown.

- Smooth KCI Drape to remove creases and seal joints ensuring an occlusive seal.
- All liners must be removed before additional drape is used.

- Press drape down around the length of the adhesive seal to ensure adhesion to the skin.

- Press down to seal all drapes.


- Minimize any creases/wrinkles in drape.
- Excess KCI Drape can be used to seal leaks.

SENSAT.R.A.C.™ PAD APPLICATION

CAUTION: Do not place the SENSAT.R.A.C.™ Pad directly over the umbilicus.

NOTE: Do not cut off the pad or insert the tubing into the foam dressing. This may occlude the tubing and cause the PREVENA PLUS™ 125 Therapy Unit to alert.


1. Choose pad application site. Give particular consideration to fluid flow and tubing positioning to allow for optional flow and avoid placement over bony prominences or within creases in the tissue.


2. Pinch drape and carefully cut an approximately **2.5 cm** hole through the drape. The hole should be large enough to allow for removal of fluid and/or exudate. It is not necessary to cut into the foam.

NOTE: Cut a hole rather than a slit, as a slit may self-seal during therapy.


3. Apply pad, which has a central disc and a surrounding outer adhesive skirt.

- Remove both backing layers 1 and 2 to expose adhesive.


- Place pad opening in central disc directly over hole in the drape.
- Apply gentle pressure on the central disc and outer skirt to ensure complete adhesion of the pad.

- Pull back on blue tab to remove pad stabilization layer.

4. Connect to the therapy unit. For connection to the PREVENA PLUS™ 125 Therapy Unit, see the **Connecting the Dressing to PREVENA PLUS™ 125 Therapy Unit** section. For connection to other KCI V.A.C.® Therapy Units, see the **Connecting the PREVEN PLUS™ Dressing with V.A.C.® Therapy Units** section.

For dressing removal, see the Dressing Removal section.

NOTE: The removable label on the supplied ruler may be used as needed to record date of dressing application or removal.


PREVENA PLUS™ CANISTER INSTALLATION

The canister used with the PREVENA PLUS™ 125 Therapy Unit is a single-use, sterile, 150 mL container with graduated markings of approximately 50 cc/mL increments.


NOTE: If the canister is not fully engaged, the PREVENA PLUS™ 125 Therapy Unit will alert.

NOTE: Only use the recommended PREVENA PLUS™ 125 Therapy Unit Canister with this product.

NOTE: Never reuse a canister.


1. Remove the canister from the sterile package.
2. Hold therapy unit and canister, vertically or horizontally, one in each hand, and slide bottom of canister into slot on bottom of therapy unit.

3. Close canister against therapy unit. The upper locking tab will click when canister is secured.


CONNECTING THE DRESSING TO THE PREVENA PLUS™ 125 THERAPY UNIT

For connecting to other KCI V.A.C.® Therapy Units, see the **Using the PREVENA PLUS™ Dressing with KCI V.A.C.® Therapy Units** section.


1. Connect the SENSAT.R.A.C.™ Pad tubing to the PREVENA PLUS™ Connector:

- Push the connectors together
- Twist the connectors to lock


2. Connect the PREVENA PLUS™ Connector to the canister by aligning and plugging connector at end of tubing onto tubing ports on side of canister. Push together firmly. Ensure clamp on tube is open. Position clamp away from patient.

3. Begin therapy.

BEGINNING THERAPY

1. Ensure the PREVENA PLUS™ Dressing has been applied as described in the **PREVENA PLUS™ CUSTOMIZABLE™ Dressing Application** section.


2. To begin therapy, press and hold center of the **On/Off** button for three seconds. The PREVENA PLUS™ 125 Therapy Unit, while in operation, may have a moderate sound emanating from the unit. All seven Therapy Life Indicators will illuminate with a green LED, indicating therapy is running.

NOTE: To interrupt therapy or turn unit off, press and hold center of the **On/Off** button for three seconds.




Once therapy is on for one hour non-stop, the seven day lifespan countdown begins and continues even if unit is turned off.

3. With therapy on, assess dressing to ensure integrity of seal.

- The dressing should have a wrinkled appearance and the foam bolster should be compressed.
- If the foam bolster is not compressed or the therapy unit alerts, see the **Alerts** section.

4. If there is any evidence of a leak, check the SENSAT.R.A.C.™ Pad and the dressing seals, tubing connectors, and canister connection, and ensure clamp is open. Refer to the **Correcting a Leak Condition** section for more information.

5. Secure excess tubing to prevent interference with patient mobility.

6. If desired, place the therapy unit into the carrying case. Ensure display is visible through the opening in the carrying case.


7. The carrying case comes with both an adjustable strap and belt clip for carrying. The belt clip and additional clips on each side and at the bottom of the carrying case provide a place where excess tubing may be wrapped and stored to help prevent/minimize tripping and strangulation.

CAUTION: Do not wear or wrap strap around neck. Do not wrap tubing around neck.

UNIT TROUBLESHOOTING

If the PREVENA PLUS™ 125 Therapy Unit will not power on, make sure batteries are charged (see **Battery Charging** section). If the therapy unit still will not turn on, contact KCI.

THERAPY LIFE INDICATORS

NOTE: Grey represents green indicators and black represents yellow indicators.

The therapy life indicators provide a visual display of the seven day therapy life cycle and the therapy life remaining. When therapy begins all seven green LEDs are illuminated. During therapy, after each 24-hour period an indicator will turn off. When eight hours of therapy time remains, the last indicator will illuminate with both a green and yellow LED simultaneously. When therapy time is about to expire, the last indicator will illuminate with a yellow LED and an alert will sound for approximately two minutes, then the therapy unit will shut off.

At the end of therapy, the therapy unit must be replaced with a new unit or alternative therapy must be used. Patients should be instructed to contact the treating physician or caregiver if therapy unit turns off and cannot be restarted before therapy is scheduled to end.

NOTE: Once therapy is on for one hour non-stop, the seven day lifespan countdown begins and continues even if the unit is turned off.

DURATION OF PREVENA PLUS™ THERAPY


- Therapy should be continuous for a minimum of two days up to a maximum of seven days.

NOTE: The PREVENA PLUS™ 125 Therapy Unit will automatically time-out after seven days of therapy.

- Patients should be instructed to contact their treating physician and not to turn therapy off unless:
 - advised by the treating physician
 - bleeding develops suddenly or in large amounts during therapy
 - there are signs of allergic reaction or infection
 - the canister is full of fluid
 - system alerts must be addressed
- Patient should be instructed to contact the treating physician if therapy unit turns off and cannot be restarted before therapy is scheduled to end, or if canister becomes full of fluid.
- At end of therapy, patient should return to treating physician for dressing removal.


ALERTS

Audible Alerts - All audible alerts will sound two beeps, escalating and repeating every 15 seconds, which will increase in volume through four cycles. The fourth cycle will produce the loudest audible beep and will repeat until the alert condition is corrected.


Alert Mute Button - Press and hold center of the **Alert Mute** button for three seconds during an alert condition to silence the audible alert for two minutes. When pressed, the **Alert Mute** button will illuminate to indicate mute has been selected. The alert will re-occur after two minutes unless the alert condition has been corrected.

The therapy unit will sound audible and display visual alerts as follows:


Alert Type	ID and Resolution
Blockage Alert 	• A solid yellow LED above the blockage symbol will turn on. • Audible blockage alert will sound two beeps repeating every 15 seconds. • When the blockage condition is resolved, audible and visual alerts will turn off. To Correct Alert Check for a full canister. Check for kinked tubing. Ensure tubing clamp is open.
Leak Alert 	• A solid yellow LED above the leak symbol will turn on. • Leak alert will sound two beeps repeating every 15 seconds. • When the leak condition is corrected, audible and visual alerts will turn off. To Correct Alert See the Correcting a Leak Condition section in this guide.
Low Battery Level Alert 	• A solid yellow LED on the battery level indicator will turn on. • Alert will sound two beeps repeating every four minutes. • A low battery alert indicates approximately two hours of therapy remain; charge batteries IMMEDIATELY to prevent disruption of therapy. • When battery is charged, audible and visual alerts will turn off. To Correct Alert Charge battery; see the Battery Charging section in this guide.
Therapy Ended 	• A solid yellow LED at the top of the Therapy Life Indicator will turn on. • The therapy unit will sound eight beeps, followed by a continuous beep for five seconds, then the therapy unit will turn off. • Notify the treating physician or caregiver: If the therapy unit has completed the seven day therapy and has timed out, and an attempt is made to turn the therapy unit on, the therapy unit will sound an alert for three seconds then shut off.
System Fault Alert  PREVENA PLUS™ 125 THERAPY UNIT	• All LEDs will turn on and flash. • Two beeps sound, repeating every 15 seconds. To Correct Alert Power the therapy unit off and then on again. If alert continues contact the treating physician or caregiver.

CORRECTING A LEAK CONDITION

When the therapy unit detects a significant leak, a visual and audible leak alert will activate (see **Alerts** section).

1. With the therapy unit on, press down through drape onto the adhesive seal all the way around to ensure seal.


2. Use fingers to smooth out any creases or wrinkles.

3. Use excess KCI Drape along the outer edge of drape to seal leaks.


CHECK CANISTER TUBING CONNECTION

1. Ensure canister is securely locked onto the therapy unit. When canister is installed, a distinct click will be heard indicating it has been properly installed.

2. Check dressing tubing connector at canister.


3. Check tubing connectors to ensure they are fully engaged and locked.

INDICATIONS THAT A LEAK CONDITION HAS BEEN CORRECTED


NOTE: Upon correcting a leak condition, a small delay will occur before the therapy unit senses the correction and silences the alerts.

The therapy unit will continue the alert until condition is corrected. When leak condition has been corrected, audible alerts will stop, and visual alerts will turn off.

The PREVENA™ Dressing will be compressed.


- Dressing compressed - system pressure acceptable.


- Dressing not compressed - system pressure not acceptable. Return to the **Correcting a Leak Condition** section to continue pressure correction steps

BATTERY CHARGING

The PREVENA PLUS™ 125 Therapy Unit is battery-operated to facilitate patient mobility. The battery charge indicator on the user interface will display three levels of charge.


- Full charge (approximately nine hours remain)
- Medium charge (approximately two - seven hours remaining)
- Low charge. When low charge is indicated approximately two hours of therapy remain. Charge unit immediately to avoid disruption of therapy

When the PREVENA PLUS™ 125 Therapy Unit is plugged into a power supply, the Power Connected icon turns yellow, indicating power is connected and system is charging. The icon will turn green when fully charged.

NOTE: Upon receipt, the PREVENA PLUS™ 125 Therapy Unit battery may not be fully charged.

NOTE: The rechargeable battery used in the PREVENA PLUS™ 125 Therapy Unit is not user accessible or replaceable.

1. Plug the mains power cord into a wall outlet.
2. Plug the other end of the mains power cord into the DC power supply.
3. Plug the DC power cord into the bottom of the therapy unit.
4. A fully discharged battery will recharge in approximately six hours.

CAUTION: Use only the charging system provided with the PREVENA PLUS™ 125 Therapy Unit. Using any other charging system may damage the therapy unit.

CAUTION: Power cords may present a tripping hazard. Ensure that power cords are out of areas where people walk.

NOTE: Power cords may have different wall plug configurations depending on country requirements.

CANISTER REMOVAL AND REPLACEMENT

1. Turn therapy off.
2. Slide dressing tubing clamp close to where tubing plugs into canister. Close clamp.
3. Unplug tubing from canister tubing ports.
4. Remove therapy unit from carrying case, if in use.
5. Press tab on canister to remove used canister from therapy unit.
6. Install new canister (see the **PREVENA PLUS™ Canister Installation** section).
7. Return therapy unit to carrying case if desired.
8. Reattach dressing tubing to canister tubing ports.
9. Release tubing clamp.
10. Turn therapy on.

NOTE: Dispose of used canister according to institution and local environmental regulations.

USING THE PREVENA PLUS™ DRESSING WITH KCI V.A.C.® THERAPY UNITS

When directed by the treating physician, PREVENA PLUS™ Dressings can be used with negative pressure wound therapy provided by ACTIV.A.C.™ and V.A.C.ULTA™ Therapy Units.


WARNING: Refer to the therapy unit's user manual for complete instructions for use and safety information before initiating therapy.

WARNING: DO NOT use with V.A.C. VERAFLU™ Therapy (Instillation) provided by the V.A.C.ULTA™ Therapy Unit. Instillation into the incision site may result in pooling of fluid which may result in maceration.

WARNING: The V.A.C.ULTA™ Therapy Unit is only indicated for use in acute care settings. Before transitioning the patient to home care, this therapy unit must be replaced with one for home use, such as the PREVENA PLUS™ 125 or ACTIV.A.C.™ Therapy Unit.

NOTE: Consider using the smallest available canister for the selected V.A.C.® Therapy Unit.

CONNECTING THE PREVENA PLUS™ DRESSING TO V.A.C.® THERAPY UNITS

1. Connect the PREVENA PLUS™ Dressing tubing to the V.A.C.® Therapy Unit canister tubing:
 - Push the connectors together
 - Twist the connectors to lock
2. Ensure clamp on canister tubing is open.

SETTING NEGATIVE PRESSURE ON THE V.A.C.® THERAPY UNITS

Set and activate V.A.C.® Therapy at -125mmHg continuous. Do not choose any other negative pressure setting or intermittent or DPC modes of negative pressure.

For V.A.C.ULTA™ Therapy Units:

- Select V.A.C.® Therapy or, when available, PREVENA™ Therapy.
- Do NOT select V.A.C.VERA FLO™ Therapy (see the **WARNING** above).
- See the **Choose Therapy** section of the V.A.C.ULTA™ User Manual for more information.

ALARM RESOLUTIONS

KCI V.A.C.® Therapy Unit alarms should be addressed in a timely manner. Refer to the appropriate therapy unit user manual for complete information on alarm resolutions. Refer to the **Correcting a Leak Condition** section for correcting a leak in the dressing.

PREVENA PLUS™ 125 THERAPY UNIT DISPOSAL

At the end of therapy, the patient should return the PREVENA PLUS™ 125 Therapy Unit to the physician for disposal. Dispose of all waste according to local requirements. Improper disposal may run the risk of regulatory non-compliance.

INSTRUCTIONS FOR PATIENTS

Review the following information with the patient prior to discharge. This information is summarized in the PREVENA PLUS™ Incision Management System Patient Guide which must be provided to the patient at discharge.

WARNING: Inform the patient that the V.A.C.ULTA™ Therapy Unit should not be sent home. They should be sent home with a KCI negative pressure therapy unit such as the PREVENA PLUS™ 125 or ACTIV.A.C.™ Therapy Units.

DAILY USE

The PREVENA PLUS™ 125 and ACTIV.A.C.™ Therapy Units are portable and small enough that they may be worn beneath clothing during normal patient activities as approved by the treating physician.

CAUTION: Advise patient to NOT SUBMERGE therapy unit or dressing in liquid and to ensure therapy unit is not pulled into a tub or sink where it may become submerged.

CAUTION: The PREVENA PLUS™ Incision Management System is a medical device not a toy. Keep away from children, pets and pests as they can damage the dressing and therapy unit and affect performance. Keep therapy unit free of dust and lint.

SLEEPING

Instruct patient to:

- place the therapy unit in a position where tubing will not become kinked or pinched.
- ensure therapy unit will not be pulled off a table or fall to the floor during sleep.

SHOWERING AND BATHING

- Do not use the PREVENA PLUS™ 125 Therapy Unit while bathing/showering or where it can fall or be pulled into a tub, shower or sink.
- Do not reach for a product that has fallen into water. Unplug the unit immediately if plugged into an electrical source. Disconnect unit from dressing and contact treating physician or caregiver.
- Light showering is permissible, bathing is not. Before showering, clamp the tubing and disconnect the dressing from the therapy unit.
- Dressing may be exposed to common shower soaps and rinsed with indirect shower stream. Do not submerge dressing. Do not remove dressing.
- When towel drying, avoid disturbing or damaging the dressing.

STRENUOUS ACTIVITY

Advise patient as to when and at what level physical activities may be resumed. It is recommended that patients avoid strenuous activity while using the PREVENA PLUS™ Incision Management System.

CLEANING

Advise patient that the therapy unit and carrying case can be wiped with a damp cloth using a mild household soap solution that does not contain bleach.

DRESSING REMOVAL

NOTE: If dressing is lifted to observe incision, do not re-adhere the same dressing; a new dressing must be applied.

WARNING: Dressings should always be removed in-line with the sutures and NEVER across the sutures.

1. Turn the therapy unit off by pressing and holding the On/Off button.

2. Gently stretch the drape/dressing horizontally to release the adhesive from the skin. Do not peel vertically. Remove the drape/dressing in-line with the sutures, NEVER across the sutures.
3. Clean any residual adhesive.

If a new dressing is to be applied:

1. Ensure the incision area is clean by using an alcohol swab or an antiseptic wipe.
2. Allow skin to dry completely.
3. Follow the **PREVENA PLUS™ CUSTOMIZABLE™ Dressing Applications** instructions.

WARNINGS AND IMPORTANT INFORMATION FOR USERS - PREVENA PLUS™ 125 THERAPY UNIT

In order for KCI products to perform properly, KCI recommends the following conditions. Failure to comply with these conditions will void any applicable warranties.

- Use this product only in accordance with these instructions and applicable product labeling.
- Assembly, operations, extensions, re-adjustments, modifications, technical maintenance or repairs must be performed by qualified personnel authorized by KCI.
- Ensure the electrical installation of the room complies with the appropriate national electrical wiring standards. To avoid the risk of electrical shock, this product must be connected to a grounded power receptacle.
- Cell phones or similar products could affect the therapy unit. Move the therapy unit away from these devices if interference is suspected.
- Do not operate this product if it has a damaged power cord, power supply or plug. If these components are worn or damaged, contact KCI.
- Do not drop or insert any object into any opening or tubing of this product.
- Do not connect this product or its components to devices not recommended by KCI.
- Do not modify the therapy unit or dressing. Do not connect the therapy unit or dressing to other devices being used.
- Use only PREVENA PLUS™ Dressings with this product.
- Keep this product away from heated surfaces.
- Equipment not suitable for use in presence of flammable anesthetic mixture with air, oxygen or nitrous oxide or in an environment in which the concentration of oxygen is: a) greater than 25% for ambient pressures up to 110 kPa; or b) the partial pressure of oxygen is greater than 27.5 kPa at ambient pressures exceeding 110 k.
- Avoid spilling fluids on any part of this product.
- Do not make any changes to the settings on the therapy unit without instructions from the treating physician.
- Small Parts - Choking Hazard
- The PREVENA PLUS™ Incision Management System is a medical device, not a toy. Keep away from children, pets and pests as they can damage the dressing and therapy unit and affect performance. Keep the therapy unit free of dust and lint as they can also damage the dressing and therapy unit and affect performance.

WARNING: The PREVENA PLUS™ 125 Therapy Unit has no serviceable parts and should not be opened, disassembled or otherwise modified by the user, and should be replaced as a unit. All assembly, operations, adjustments, modifications, maintenance and repairs must be carried out by qualified personnel authorized by KCI.

Electric Shock Hazard - Do not open any electrical cover on the therapy unit. There are no serviceable parts. Refer to qualified KCI service personnel.

Fluids remaining on the electronic controls can cause corrosion that may cause the electronic components to fail. Component failures may cause the unit to operate erratically, possibly producing potential hazards to patient and staff. If spills do occur, unplug the unit immediately and clean with an absorbent cloth. Ensure there is no moisture in or near the power connection and power supply components before reconnecting power. If the product does not work properly, contact KCI.

CUSTOMER CONTACT INFORMATION

For questions regarding this product, supplies, maintenance or additional information about KCI products and services, please contact KCI or a KCI authorized representative, or:

In the US call 1-800-275-4524 or visit www.acelity.com

Outside the US visit www.kci-medical.com

KCI USA, Inc. 12930 IH 10 West, San Antonio, TX 78249

PREVENA PLUS™ THERAPY SYSTEM SPECIFICATIONS

Dimensions 60300EP:	3.5" W x 6.4" H x 2.16" D (8.9 x 16.3 x 5.49 cm)
Dimensions 60300P:	3.5" W x 6.4" H x 1.81" D (8.9 x 16.3 x 4.62 cm)
Weight 60300EP (with empty canister attached):	~.64 lbs (~.29 kg)
Weight 60300P (with empty canister attached):	~.75 lbs (~.34 kg)
Pressure:	125 mmHg (16.7 kPa)
Canister Volume:	150 mL

Electrical:

Battery Run Life:	~8.5 hours
Battery Charge Time:	~6 hours from a fully discharged state
External Power Supply Input:	100-240VAC 0.5A-0.3A 50 - 60 Hz
External Power Supply Output 60300EP:	5V, 2.0 A
External Power Supply Output 60300P:	9V, 2.2 A
Patient and Enclosure Leakage Current:	<100 Microamps

Environmental Conditions:

Storage/Transport Conditions

Temperature Range:	0°F (-18°C) to 140°F (60°C)
Relative Humidity Range:	15-95% non-condensing

Operating Conditions

Temperature Range:	41°F (5°C) to 104°F (40°C)
Relative Humidity Range:	15-93% non-condensing
Atmospheric Pressure:	1060 hpa (-1253 ft/-381. 9 m) to 700 hpa (9878 ft/3010 m)
Expected Service Life	7.5 days

IEC Classification

Medical Equipment

Type BF, Applied Part

Class II

IP22 - Protection against solid objects greater than 12.5 mm and against liquid water falling for short periods of time.

The dressing components of the PREVENA PLUS™ 125 Incision Management System are considered Applied Parts under IEC 60601-1.

Conforms to AAMI ES60601-1:2005 +A1:2012, IEC 60601-1-6 Ed.3.1:2013, IEC 60601-1-8 Ed.2.1:2012, IEC 60601-1-11 Ed.2.0:2015, CSA C22.2#60601-1:2014 Ed.3

PREVENA PLUS™ 125 THERAPY UNIT ELECTROMAGNETIC COMPATIBILITY

The following are guidance and manufacturer's declarations regarding EMC for the PREVENA PLUS™ 125 Therapy Unit.

- The PREVENA PLUS™ 125 Therapy Unit needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the following pages.

WARNING: This equipment is intended for use by healthcare professionals only. As with all electrical medical equipment, this equipment may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures such as re-orienting or relocating the PREVENA PLUS™ 125 Therapy Unit or shielding the location.

- Portable and Mobile RF communications equipment, RFID readers, electronic article surveillance (anti-theft) equipment and metal detectors can affect the performance of the PREVENA PLUS™ 125 Therapy Unit. Please use the guidelines and recommendations specified in Tables 204 and 206.
- Other medical equipment or systems can produce electromagnetic emissions and therefore can interfere with the functionality of the PREVENA PLUS™ 125 Therapy Unit. Care should be used when operating the PREVENA PLUS™ 125 Therapy Unit adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the PREVENA PLUS™ 125 Therapy Unit and the other equipment should initially be observed to verify normal operation in the configuration in which it will be used.
- The electrical cables, external power supplies and accessories listed or referenced in this manual have been shown to comply with the test requirements listed in the following tables. Care should be taken to use only manufacturer-recommended cables, power supplies and accessories with the PREVENA PLUS™ 125 Therapy Unit. If a third-party supplier offers cables, external power supplies and electrical accessories for use with the PREVENA PLUS™ 125 Therapy Unit and they are not listed or referenced in this manual, it is the responsibility of that third-party supplier to determine compliance with the standards and tests in the following tables.
- The use of electrical cables and accessories other than those specified in this manual or referenced documents may result in increased electromagnetic emissions from the PREVENA PLUS™ 125 Therapy Unit or decreased electromagnetic immunity of the PREVENA PLUS™ 125 Therapy Unit.
- Portable and mobile RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the PREVENA PLUS™ 125 Therapy Unit including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. More precisely, the minimum recommended separation distance should be calculated from the equation applicable to the frequency of the transmitter, as noted in the guidance below.
- NOTE: This equipment has been tested and found to comply with the limits for medical devices to IEC 60601-1-2: 2014 4th edition. These limits and test levels are intended to provide reasonable safety with regard to electromagnetic disturbances when the device is used in a typical medical installation.

Table 201

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The PREVENA PLUS™ 125 Therapy Unit is intended for use in the electromagnetic environment specified below. The customer or the end user of the PREVENA PLUS™ 125 Therapy Unit should assure that it is used in such an environment.		
Emission Test	Compliance	Electromagnetic Environment - Guidance
RF emissions - CISPR 11 (Radiated & Conducted)	Group 1	The PREVENA PLUS™ 125 Therapy Unit uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The PREVENA PLUS™ 125 Therapy Unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions - CISPR 11 (Radiated & Conducted)	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Table 202

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The PREVENA PLUS™ 125 Therapy Unit is intended for use in the electromagnetic environment specified below. The customer or the end user of the unit should assure it is used only in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Intended Electromagnetic Environment
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines 100 kHz repetition frequency ±1 kV for input/output lines 100 kHz repetition frequency	±2 kV for power supply lines 100 kHz repetition frequency ±1 kV for input/output lines 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or home healthcare environment.
Surge IEC 61000-4-5	±1 kV differential mode (line - line) ±2 kV common mode (line - earth)	±1 kV differential mode (line - line) ±2 kV common mode (line - earth)	Mains power quality should be that of a typical commercial or home healthcare environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Dips: 0% Ut for 1 cycle 70% Ut for 25 cycles at 50 Hz or 30 cycles at 60 Hz Single phase: at 0° Interruptions: 0% Ut for 250 cycles at 50 Hz or 300 cycles at 60 Hz	Dips: 0% Ut for 1 cycle 70% Ut for 25 cycles at 50 Hz or 30 cycles at 60 Hz Single phase: at 0° Interruptions: 0% Ut for 250 cycles at 50 Hz or 300 cycles at 60 Hz	Product has internal battery backup. If the user of the PREVENA PLUS™ 125 Therapy Unit requires continued operation during power mains interruptions, it is recommended that the unit be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or home healthcare environment.
NOTE: Ut is the A.C. mains voltage prior to application of the test level.			

Table 204

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The PREVENA PLUS™ 125 Therapy Unit is intended for use in the electromagnetic environment specified below. The customer or the end user of the unit should assure it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3Vrms 150kHz to 80 MHz 6Vrms in ISM and amateur radio bands between 150 kHz and 80 MHz 80% AM at 1kHz 10 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz	3Vrms 150kHz to 80 MHz 6Vrms in ISM and amateur radio bands between 150 kHz and 80 MHz 80% AM at 1kHz 10 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the PREVENA PLUS™ 125 Therapy Unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 0.35\sqrt{P}$ 80MHz to 800 MHz $d = 0.7\sqrt{P}$ 800MHz to 2.7GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended minimum separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1: At 80MHz and 800MHz, the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from objects, structures and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the PREVENA PLUS™ 125 Therapy Unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the PREVENA PLUS™ 125 Therapy Unit.			

Table 205

Guidance and Manufacturer's Declaration - Electromagnetic Immunity						
Test Frequency (MHz)	Band ^a (MHz)	Service ^a	Modulation ^b	Maximum Power (W)	Distance (m)	Immunity Test Level (V/m)
385	380 - 390	TETRA 400	Pulse modulation ^b 18 Hz	1,8	0,3	27
450	430 - 470	GMRS 460, FRS 460	FM ^c ±5 kHz deviation 1 kHz sine	2	0,3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation ^b 217 Hz	0,2	0,3	9
745						
780						
810	800 - 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^b 18 Hz	2	0,3	28
870						
930						
1 720	1 700 - 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3 4, 25; UMTS	Pulse modulation ^b 217 Hz	2	0,3	28
1 845						
1 970						
2 450	2 400 - 2 570	Bluetooth WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^b 217 Hz	2	0,3	28
5 240	5 100 - 5 800	WLAN 802.11 a/n	Pulse modulation ^b 217 Hz	0,2	0,3	9
5 500						
5 785						
NOTE: If necessary to achieve the Immunity Test Level, the distance between the transmitting antenna and the ME Equipment or ME System may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.						
^a For some services, only the uplink frequencies are included.						
^b The carrier shall be modulated using a 50% duty cycle square wave signal.						
^c As an alternative to FM modulation, a 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.						

Table 206

Recommended separation distances between portable and mobile RF communications equipment and the PREVENA PLUS™ 125 Therapy Unit			
The PREVENA PLUS™ 125 Therapy Unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the PREVENA PLUS™ 125 Therapy Unit can help prevent electromagnetic interference by maintaining a minimum distance between the portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter in watts (W)	Separation distance according to frequency of transmitter in meters (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0.01	.12	.04	.07
0.1	.38	.11	.22
1	1.2	.35	.7
10	3.8	1.1	2.2
100	12	3.5	7
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

INCLUDED POWER SUPPLIES

Part Number	Description	Manufacturer	Max Length
44000919	Power Supply, 240V	Inventus Power	3.10 m
44001674	Power Supply, 5V, 10W	Inventus Power	1.03 m
413628	Cord, V.A.C.* Power, US	Consolidated Wire	2.08 m
413625	Cord, V.A.C.* Power, UK-240V	Consolidated Wire	2.08 m
413626	Cord, V.A.C.* Power, EU-240V	Consolidated Wire	2.08 m
413627	Cord, V.A.C.* Power, AU/NZ 240V	Consolidated Wire	2.08 m
413992	Cord, V.A.C.* Power, CH-240V	Consolidated Wire	2.08 m
414961	Cord, V.A.C.* Power, South Africa	Consolidated Wire	2.08 m

The use of electrical cables and accessories other than those specified in the supplied instructions for use or referenced documents may result in increased electromagnetic emissions from the PREVENA PLUS™ 125 Therapy Unit or decreased electromagnetic immunity of the PREVENA PLUS™ 125 Therapy Unit.

BIBLIOGRAPHY OF PUBLISHED STUDIES

- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 2009 June 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 Feb;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 2013 December 1;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 2013 January 30.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 2013 Feb 19. [Epub ahead of print]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 2013 December 1;56(12):1403-8.
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 2011 December 6;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery 2013 Oct 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 2012 Aug 21(8):386-8.
- Faruqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 2013 September 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 2012 Nov 24(11):308-11.
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 Nov-Dec;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. 2013 May;145(5):1387-92.
- Haghsheenas Kashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. 2011 Jan/Feb;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 2013 September 20;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 2013 December 1;2(12):276-84.

- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. 2011 Sep;19(5):588-96.
- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations. 2013 September 20.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. 2012 Jun;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. 2013 Mar;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. 2012 Apr;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. 2013 Jun;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. 2010 Jun;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. 2006 Jun;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. 2012 Jan;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. 2013 Mar;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. 2012 Dec;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. 2012 Mar;19(1):67-75.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 2009 juni 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 feb;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 2013 december 1;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 2013 januari 30.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 2013 feb 19. [Epub voorafgaand aan gedrukte uitgave]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 2013 december 1;56(12):1403-8.
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 2011 december 6;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery 2013 oktober 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 2012 augustus 21(8):386-8.
- Faruqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 2013 september 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 2012 november 24(11):308-11.
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 nov-dec;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. 2013 mei;145(5):1387-92.
- Haghsheenas Kashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. 2011 jan/feb;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 2013 september 20;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 2013 december 1;2(12):276-84.

- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. 2011 september;19(5):588-96.
- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations. 2013 september 20.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. 2012 juni;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. 2013 maart;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. 2012 april;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. 2013 juni;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. 2010 juni;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. 2006 juni;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. 2012 januari;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. 2013 maart;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. 2012 dec;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. 2012 maart;19(1):67-75.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. Juin 2009 ; 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. Fév. 2011 ; 8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 1 décembre 2013 ; 12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 30 janvier 2013.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 19 fév. 2013 [Publication électronique avant impression]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 1 décembre 2013 ; 56(12):1403-8.
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 6 décembre 2011 ; 6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery. Oct. 2013 ; 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 21 août 2012 ; (8):386-8.
- Faruqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 1 septembre 2013 ; 5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 24 nov 2012 ; (11):308-11.
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. Nov./Déc. 2006 ; 20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. Mai 2013 ; 145(5):1387-92.
- Haghsheenasakshani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. Jan./Fév. 2011 ; 11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 2013 September 20;13:e49.20 septembre 2013 ;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 1 décembre 2013 ; 2(12):276-84.

- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. Sep. 2011 ; 19(5):588-96.
- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations. 2013 September 20.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. Jun 2012 ; 255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. Mar. 2013 ; 57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. Avr. 2012 ; 36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. Jun 2013 ; 14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. Jun 2010 ; 19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. Jun 2006 ; 60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. Jan. 2012 ; 26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. Mar 2013 ; 66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. Déc. 2012 ; 204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. 2012 Mar;19(1):67-75.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 2009 June 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 Feb;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 2013 December 1;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 2013 January 30.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 2013 Feb 19. [Epub ahead of print]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 2013 December 1;56(12):1403-8.
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 2011 December 6;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery. 2013 Oct 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 2012 Aug 21(8):386-8.
- Farooqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 2013 September 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 2012 Nov 24(11):308-11.
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 Nov-Dec;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. 2013 May;145(5):1387-92.
- Haghshenas Kashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. 2011 Jan/Feb;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty. 2013 September 20;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 2013 December 1;2(12):276-84.
- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. 2011 Sep;19(5):588-96.

- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations. 2013 September 20.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. 2012 Jun;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. 2013 Mar;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. 2012 Apr;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. 2013 Jun;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. 2010 Jun;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. 2006 Jun;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. 2012 Jan;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. 2013 Mar;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. 2012 Dec;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. 2012 Mar;19(1):67-75.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. Junio de 2009; 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy". International Wound Journal. Febrero de 2011; 8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures". Techniques in Foot & Ankle Surgery. 1 de diciembre de 2013; 12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy". American Journal of Surgery. 30 de enero de 2013.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report". International Wound Journal. 19 de febrero de 2013 [publicación electrónica antes de su impresión]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery". Diseases of the Colon & Rectum. 1 de diciembre de 2013; 56(12):1403-8.
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients". Journal of Cardiothoracic Surgery. 6 de diciembre de 2011; 6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study". Annals of Plastic Surgery. Octubre de 2013 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown". Journal of Wound Care. Agosto de 2012; 21(8):386-8.
- Faruqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively". Journal of Diabetic Foot Complications. 1 de septiembre de 2013; 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study". Wounds. Noviembre de 2012; 24(11):308-11.
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy". Journal of Orthopaedic Trauma. Nov-dic de 2006; 20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy". Journal of Thoracic and Cardiovascular Surgery. Mayo de 2013; 145(5):1387-92.
- Haghsheenas Kashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass". The British Journal of Diabetes and Vascular Disease. Ene-feb de 2011; 11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review". Eplasty. 20 de septiembre de 2013; 13:e49.
- Karlakki S. et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication". Bone and Joint Research. 1 de diciembre de 2013; 2(12):276-84.

- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system". Wound Repair and Regeneration. Septiembre de 2011;19(5):588-96.
- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study". Surgical Innovations 20 de septiembre de 2013.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study". Annals of Surgery. Junio de 2012;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients". Journal of Vascular Surgery. Marzo de 2013; 57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty". International Orthopaedics. Abril de 2012;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions". Surgical Infections. Junio de 2013;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery". Journal of Surgical Orthopaedic Advances. Junio de 2010; 19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma". Journal of Trauma. Junio de 2006; 60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures". Journal of Orthopedic Trauma. Junio de 2012; 26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications". Journal of Plastic, Reconstructive, and Aesthetic Surgery. Marzo de 2013; 66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures". American Journal of Surgery. Diciembre de 2012; 204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics". Surgical Innovations. Marzo de 2012; 19(1):67-75.

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- Atkins BZ, v.d. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 16 Haziran 2009(2):140-6.
- Atkins BZ, v.d. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 Şubat;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 1 Aralık 2013;12(4):172-4.
- Blackham AU, v.d. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 30 Ocak 2013.
- Bollero D, v.d. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 19 Şubat 2013. [Epub baskı öncesi]
- Bonds AM, v.d. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 1 Aralık 2013;56(12):1403-8
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 6 Aralık 2011;6(1):160.
- Conde-Green A, v.d. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery 2013 Ekim 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 2012 Ağustos 21(8):386-8.
- Faruqi L, v.d. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 1 Eylül 2013;5(2):44-7.
- Glaser DA, v.d. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 2012 Kasım 24(11):308-11
- Gomoll AH, v.d. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 Kasım-Aralık;20(10):705-9.
- Grauhan O, v.d. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. 2013 Mayıs;145(5):1387-92.
- Haghshenas-kashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. 2011 Ocak/Şubat;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 20 Eylül 2013;13:e49.
- Karlakki S, v.d. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 1 Aralık 2013;2(12):276-84.
- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. 2011 Eylül;19(5):588-96.

- Mark KS, v.d. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A Pilot Study." Surgical Innovations 20 Eylül 2013.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. 2012 Haziran;255(6):1043-7.
- Matatov T, v.d. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. 2013 Mart;57(3):791-5.
- Pachowsky M, v.d. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. 2012 Nisan;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. 2013 Haziran;14(3):270-4.
- Reddix RN, v.d. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. 2010 Haziran;19(2):91-7.
- Stannard JP, v.d. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. 2006 Haziran;60(6):1301-6.
- Stannard JP, v.d. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. 2012 Ocak;26(1):37-42.
- Tauber R, v.d. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. 2013 Mart;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. 2012 Aralık;204(6):1021-4.
- Wilkes RP, v.d. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. 2012 Mart;19(1):67-75.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 2009 June 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 Feb;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 2013 December 1;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 2013 January 30.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 2013 Feb 19. [Εpub πριν από την εκτύπωση]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 2013 December 1;56(12):1403-8.
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 2011 December 6;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery 2013 Oct 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 2012 Aug 21(8):386-8.
- Faruqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 2013 September 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 2012 Nov 24(11):308-11.
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 Nov-Dec;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. 2013 May;145(5):1387-92.
- Haghsheenas Kashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. 2011 Jan/Feb;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 2013 September 20;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 2013 December 1;2(12):276-84.

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- Atkins BZ, et al. "Does Negative Pressure Wound Therapy Have a Role in Preventing Poststernotomy Wound Complications?" Surgical Innovations. 2009 June 16(2):140-6.
- Atkins BZ, et al. "Laser Doppler flowmetry assessment of peristernal perfusion after cardiac surgery: beneficial effect of negative pressure therapy." International Wound Journal. 2011 Feb;8(1):56-62.
- Berkowitz MJ. "Use of a Negative Pressure Incisional Dressing After Surgical Treatment of Calcaneal Fractures." Techniques in Foot & Ankle Surgery. 2013 December 1;12(4):172-4.
- Blackham AU, et al. "Prevention of surgical site infections in high-risk patients with laparotomy incisions using negative-pressure therapy." American Journal of Surgery. 2013 January 30.
- Bollero D, et al. "Negative pressure surgical management after pathological scar surgical excision: a first report." International Wound Journal. 2013 Feb 19. [sähköinen julkaisu ennen painettua lehteä]
- Bonds AM, et al. "Incisional Negative Pressure Wound Therapy Significantly Reduces Surgical Site Infection in Open Colorectal Surgery." Diseases of the Colon & Rectum. 2013 December 1;56(12):1403-8.
- Colli A. "First Experience With a New Negative Pressure Incision Management System on Surgical Incisions After Cardiac Surgery in High Risk Patients." Journal of Cardiothoracic Surgery. 2011 December 6;6(1):160.
- Conde-Green A, et al. "Incisional Negative-Pressure Wound Therapy versus conventional dressings following abdominal wall reconstruction. A comparative study." Annals of Plastic Surgery 2013 Oct 71(4):394-7.
- Dutton M, Curtis K. "Well-wound therapy: use of NPWT to prevent laparotomy breakdown." Journal of Wound Care. 2012 Aug 21(8):386-8.
- Faruqi L, et al. "Use of an Incision-Line Negative Pressure Wound Therapy Technique to Protect High-Risk Diabetic Foot Wounds, Postoperatively." Journal of Diabetic Foot Complications. 2013 September 1;5(2):44-7.
- Glaser DA, et al. "Negative pressure therapy for closed spine incisions: A pilot study." Wounds. 2012 Nov 24(11):308-11.
- Gomoll AH, et al. "Incisional vacuum-assisted closure therapy." Journal of Orthopaedic Trauma. 2006 Nov-Dec;20(10):705-9.
- Grauhan O, et al. "Prevention of poststernotomy wound infections in obese patients by negative pressure wound therapy." Journal of Thoracic and Cardiovascular Surgery. 2013 May;145(5):1387-92.
- Haghsheenas Kashani A, Varcoe RL. "A new negative pressure dressing (Prevena) to prevent wound complications following lower limb distal arterial bypass." The British Journal of Diabetes and Vascular Disease. 2011 Jan/Feb;11(1):21-4.
- Ingargiola MJ, Daniali LN, Lee ES. "Does the application of incisional negative pressure therapy to high-risk wounds prevent surgical site complications? A systematic review." Eplasty 2013 September 20;13:e49.
- Karlakki S, et al. "Negative pressure wound therapy for management of the surgical incision in orthopaedic surgery: A review of evidence and mechanisms for an emerging indication." Bone and Joint Research. 2013 December 1;2(12):276-84.

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- Kilpadi DV, Cunningham MR. "Evaluation of closed incision management with negative pressure wound therapy (CIM): hematoma/seroma and involvement of the lymphatic system." Wound Repair and Regeneration. 2011 Sep;19(5):588-96.
- Mark KS, et al. "Incisional Negative Pressure Therapy to Prevent Wound Complications Following Cesarean Section in Morbidly Obese Women: A pilot study." Surgical Innovations. 2013 September 20.
- Masden D, Goldstein J, Endara M, Xu K, Steinberg J, Attinger C. "Negative pressure wound therapy for at-risk surgical closures in patients with multiple comorbidities: a prospective randomized controlled study." Annals of Surgery. 2012 Jun;255(6):1043-7.
- Matatov T, et al. "Experience with a new negative pressure incision management system in prevention of groin wound infection in vascular surgery patients." Journal of Vascular Surgery. 2013 Mar;57(3):791-5.
- Pachowsky M, et al. "Negative pressure wound therapy to prevent seromas and treat surgical incisions after total hip arthroplasty." International Orthopaedics. 2012 Apr;36(4):719-22.
- Pauli EM, Krpata DM, Novitsky YW, Rosen MJ. "Negative Pressure Therapy for High-Risk Abdominal Wall Reconstruction Incisions." Surgical Infections. 2013 Jun;14(3):270-4.
- Reddix RN, et al. "The effect of incisional negative pressure therapy on wound complications after acetabular fracture surgery." Journal of Surgical Orthopaedic Advances. 2010 Jun;19(2):91-7.
- Stannard JP, et al. "Negative pressure wound therapy to treat hematomas and surgical incisions following high-energy trauma." Journal of Trauma. 2006 Jun;60(6):1301-6.
- Stannard JP, et al. "Incisional negative pressure wound therapy after high-risk lower extremity fractures." Journal of Orthopedic Trauma. 2012 Jan;26(1):37-42.
- Tauber R, et al. "Inguinal lymph node dissection: epidermal vacuum therapy for prevention of wound complications." Journal of Plastic, Reconstructive, and Aesthetic Surgery. 2013 Mar;66(3):390-6.
- Vargo D. "Negative pressure wound therapy in the prevention of wound infection in high risk abdominal wound closures." American Journal of Surgery. 2012 Dec;204(6):1021-4.
- Wilkes RP, et al. "Closed incision management with negative pressure wound therapy (CIM): biomechanics." Surgical Innovations. 2012 Mar;19(1):67-75.

EN - SYMBOLS USED, DE - VERWENDETE SYMBOLE, NL - GEBRUIKTE SYMBOLEN, FR - SYMBOLES UTILISÉS, IT - SIMBOLI UTILIZZATI, ES - SÍMBOLOS UTILIZADOS, DA - ANVENDETE SYMBOLER, SV - SYMBOLER SOM ANVÄNDS, PTBR - SÍMBOLOS USADOS, TR - KULLANILAN SEMBOLLER, EL - ΣΥΜΒΟΛΑ ΠΟΥ ΧΡΗΣΙΜΟΠΟΙΟΥΝΤΑΙ, FI - KÄYTETYT SYMBOLIT, NO - SYMBOLER SOM BRUKES

EN - Refer to Clinician Guide
DE - Siehe „Handbuch für den Anwender“
NL - Raadpleeg de handleiding voor behandelaars
FR - Consulter le Guide du clinicien
IT - Fare riferimento alla Guida per il medico
ES - Consultar la Guía de referencia para facultativos
DA - Se den kliniske vejledning
SV - Se den kliniska guiden
PTBR - Consultar o guia da clínica
TR - Klinisyen Kılavuzuna Bakın
EL - Ανατρέξτε στον Κλινικό οδηγό
FI - Lue käyttöopas hoitavalle lääkärille
NO - Se den medisinske veiledningen

EN - Consult Instructions for Use
DE - Gebrauchsanweisung beachten.
NL - Raadpleeg de gebruiksaanwijzing
FR - Consulter le mode d'emploi
IT - Consultare le Istruzioni per l'uso
ES - Consultar las Instrucciones de uso
DA - Læs brugervejledningen
SV - Se bruksanvisning
PTBR - Consulte as instruções de uso
TR - Kullanım Talimatlarına Başvurun
EL - Συμβουλευθείτε τις οδηγίες χρήσης
FI - Lue käyttöohjeet
NO - Se bruksanvisningen

LOT

EN - Lot Number
DE - Chargennummer
NL - Partijnummer
FR - Numéro du lot
IT - Numero di lotto
ES - Número de lote
DA - Lotnummer
SV - Lotnummer
PTBR - Número do lote
TR - Lot Numarası
EL - Αριθμός παρτίδας
FI - Eränumero
NO - Lotnummer

EN - Type BF applied part
DE - Anwendungsteil Typ BF
NL - Type BF toegepast onderdeel
FR - Partie appliquée de type BF
IT - Parte applicata di tipo BF
ES - Pieza aplicada de tipo BF
DA - Anvendt del type BF
SV - Applicerad del av typ BF
PTBR - Peça aplicada tipo BF
TR - BF Tipi uygulamlı parça
EL - Εφαρμοζόμενο μέρος τύπου BF
FI - Potilasta koskettava BF-tyyppin osa
NO - Anvendt del av typen BF

STERILE R

EN - Sterile using radiation
DE - Steril durch Bestrahlung
NL - Gesteriliseerd door middel van straling
FR - Stérilisation par irradiation
IT - Sterilizzato con radiazioni
ES - Esterilización por radiación
DA - Strålesteriliseret
SV - Strålsteriliserad
PTBR - Esterilizada por Radiação
TR - Radyasyon sterilizasyonu
EL - Αποστειρωμένο με χρήση ακτινοβολίας
FI - Steriloitu säteilyttämällä
NO - Steril ved hjelp av stråling

IP22

EN - Ingress Protection
DE - Schutz gegen Eindringen
NL - Bescherming tegen binnendringen
FR - Niveau de protection
IT - Protezione d'ingresso
ES - Protección contra penetración
DA - Kapslingsklasse
SV - Skyddskapsling
PTBR - Proteção de entrada
TR - Ingress Protection
EL - Βαθμός προστασίας IP
FI - Suojaustaso
NO - Beskyttelse mot inntrengning

EN - Date of Manufacture
DE - Herstellungsdatum
NL - Productiedatum
FR - Date de fabrication
IT - Data di produzione
ES - Fecha de fabricación
DA - Fremstillingsdato
SV - Tillverkningsdatum
PTBR - Data de fabricação
TR - Üretim Tarihi
EL - Ημερομηνία κατασκευής
FI - Valmistuspäivä
NO - Produksjonsdato

EN - Manufacturer
DE - Hersteller
NL - Fabrikant
FR - Fabricant
IT - Produttore
ES - Fabricante
DA - Producent
SV - Tillverkare
PTBR - Fabricante
TR - Üretici
EL - Κατασκευαστής
FI - Valmistaja
NO - Produsent

EN - Content Information
DE - Inhaltsinformationen
NL - Gegevens over de inhoud
FR - Contenu de l'emballage
IT - Informazioni sul contenuto
ES - Información de contenido
DA - Kontaktplysninger
SV - Innehållsinformation
PTBR - Informações do conteúdo
TR - İçerik Bilgisi
EL - Πληροφορίες περιεχομένου
FI - Tietoa sisällöstä
NO - Informasjon om innhold

EN - Use By
DE - Verfallsdatum
NL - Uiterste gebruiksdatum
FR - Date limite d'utilisation
IT - Data di scadenza
ES - Fecha de caducidad
DA - Anvendes før
SV - Använd före
PTBR - Usar até
TR - Son Kullanma Tarihi
EL - Ημερομηνία λήξης
FI - Viimeinen käyttöpäivä
NO - Brukes innen

EN - Fragile
DE - Zerbrechlich
NL - Breekbaar
FR - Fragile
IT - Fragile
ES - Frágil
DA - Forsigtig
SV - Ömtålig
PTBR - Frágil
TR - Kırılabilir Eşya
EL - Εύθραστο
FI - Särkyvä
NO - Skjør

EN - Catalog Number
DE - Katalog-Nummer
NL - Catalogusnummer
FR - Référence
IT - Numero di catalogo
ES - Número de catálogo
DA - Katalognummer
SV - Katalognummer
PTBR - Número de catálogo
TR - Katalog Numarası
EL - Αριθμός καταλόγου
FI - Tuotenumero
NO - Katalognummer

EN - Keep Dry
DE - Trocken lagern.
NL - Droog houden
FR - Protéger de l'humidité
IT - Conservare in un luogo asciutto
ES - Mantener seco
DA - Opbevares tørt
SV - Förvaras torrt
PTBR - Proteger contra umidade
TR - Kuru Tutun
EL - Διατηρείτε το προϊόν στεγνό
FI - Pidettävä kuivana
NO - Må holdes tørt

EN - Single Use Only
DE - Nur für den Einmalgebrauch
NL - Uitsluitend voor eenmalig gebruik
FR - Usage unique
IT - Esclusivamente monouso
ES - Un solo uso
DA - Kun til engangsbrug
SV - Endast för engångsbruk
PTBR - Usar apenas uma vez
TR - Yalnızca Tek Kullanımlıdır
EL - Για μία μόνο χρήση
FI - Kertakäyttöinen
NO - Kun til engangsbruk

EN - MR Unsafe
DE - Nicht für die Verwendung im MRT geeignet
NL - MR-onveilig
FR - Non protégée contre le champ magnétique de la RM
IT - Non sicuro per la risonanza magnetica
ES - No es seguro en un entorno magnético
DA - Ikke MR-sikker
SV - Inte MR-säker
PTBR - Não é seguro sob RM
TR - MR Ortamında Kullanımı Güvenli Değildir
EL - Μη ασφαλές σε περιβάλλον διεξαγωγής μαγνητικής τομογραφίας (MRI)
FI - Ei sovellu magneettikuvaukseen (MR Unsafe)
NO - Ikke MR-sikkert

EN - Tripping Hazard
DE - Stolpergefahr
NL - Struikelgevaar
FR - Risque de chute
IT - Rischio di inciampo
ES - Riesgo de tropiezo
DA - Snubelfare
SV - Snubbelrisk
PTBR - Risco de tropeço
TR - Takılma Tehlikesi
EL - Κίνδυνος παραπατήματος
FI - Kompastumisvaara
NO - Snubelfare

EN - Do Not Resterilize
DE - Nicht wieder sterilisieren.
NL - Niet opnieuw steriliseren
FR - Ne pas stériliser une seconde fois
IT - Non risterilizzare
ES - No reesterilizar
DA - Må ikke steriliseres igen
SV - Omsteriliseras ej
PTBR - Não reesterilizar
TR - Yeniden Sterilize Etmeyin
EL - Μην επαναποστειρώνετε
FI - Älä steriloi uudelleen
NO - Ikke steriliser på nytt

EN - Class II Device
DE - Gerät der Schutzklasse II
NL - Klasse II-apparaat
FR - Équipement de classe II
IT - Dispositivo di classe II
ES - Dispositivo de Clase II
DA - Klasse II-udstyr
SV - Utrustning av klass II
PTBR - Dispositivo Classe II
TR - Sınıf II Cihaz
EL - Συσκευή κατηγορίας II
FI - Luokan II laite
NO - Utstyr i klasse II

EN - No Bathing or Showering
DE - Nicht Baden oder Duschen
NL - Niet baden of douchen
FR - Bain et douche interdits
IT - Divieto di utilizzo durante il bagno o la doccia
ES - No bañarse ni ducharse
DA - Må ikke tages med i karbad eller brusebad
SV - Ej bad eller dusch
PTBR - Não lavar em chuveiro ou imersão
TR - Banyo Yapmayın veya Duş Almayın
EL - Όχι χρήση στο μπάνιο ή το ντους
FI - Älä käytä kylvyn tai suihkun aikana
NO - Bading eller dusjing forbudt

EN - Temperature Limits
DE - Temperaturgrenzen
NL - Temperatuurgrenzen
FR - Limites de température
IT - Limiti di temperatura
ES - Límites de temperatura
DA - Temperaturområde
SV - Temperaturgränser
PTBR - Limites de temperatura
TR - Sıcaklık Sınırları
EL - Όρια θερμοκρασίας
FI - Lämpötilarajat
NO - Temperaturgrenser

EN - Do not use if package is damaged or open
DE - Nicht verwenden, wenn die Verpackung beschädigt ist oder bereits geöffnet wurde.
NL - Niet gebruiken als de verpakking is beschadigd of geopend
FR - Ne pas utiliser si l'emballage est endommagé ou ouvert
IT - Non utilizzare se la confezione è danneggiata o aperta
ES - No usar si el embalaje estuviera dañado o roto
DA - Må ikke anvendes, hvis emballagen er beskadiget eller har været åbnet
SV - Innehållsinformation
PTBR - Não usar se a embalagem estiver danificada ou aberta
TR - Hasar görmüş veya açılmış ambalajı kullanmayın
EL - Μην το χρησιμοποιείτε εάν η συσκευασία έχει ανοιχτεί ή υποστεί ζημιά
FI - Tuotetta ei saa käyttää, jos pakkaus on vaurioitunut tai auki
NO - Skal ikke brukes hvis pakken er skadet eller åpnet

EN - Contains Phthalates (SENSAT.R.A.C.TM Pad Tubing)
DE - Enthält Phthalate (SENSAT.R.A.C.TM-Pad-Schlauch)
NL - Bevat ftalaten (slang van de SENSAT.R.A.C.TM Pad)
FR - Contient des phtalates (tubulure du tampon SENSAT.R.A.C.TM)
IT - Contiene ftalati (tubo del pad SENSAT.R.A.C.TM)
ES - Contiene ftalatos (tubo de la Interfase SENSAT.R.A.C.TM)
DA - Indeholder phthalater (slangen til SENSAT.R.A.C.TM Pad)
SV - Innehåller ftalater (SENSAT.R.A.C.TM-padslang)
PTBR - Contém Ftalatos (Tubulação do Coletor SENSAT.R.A.C.TM)
TR - Ftalat içerir (SENSAT.R.A.C.TM Ped Hortumu)
EL - Περιέχει φθαλικές ενώσεις (σωλήνωση επιθέματος SENSAT.R.A.C.TM)
FI - Sisältää ftalaatteja (SENSAT.R.A.C.TM-tyynyn letku)
NO - Inneholder ftalater (SENSAT.R.A.C.TM-puteslange)

EN - Authorized Representative in the European Community
DE - Autorisierte Vertretung in der EU
NL - Geautoriseerde vertegenwoordiger in de Europese Gemeenschap
FR - Mandataire européen
IT - Rappresentante autorizzato per l'Unione Europea
ES - Representante autorizado en la Comunidad Europea
DA - Godkendt repræsentant i EU
SV - Auktoriserad representant i Europeiska unionen
PTBR - Representante Autorizado na Comunidade Europeia
TR - Avrupa Topluluğundaki Yetkili Temsilci
EL - Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
FI - Valtuutettu edustaja Euroopan yhteisön alueella
NO - Autorisert representant i EU

EN - Caution
DE - Achtung
NL - Let op
FR - Attention
IT - Attenzione
ES - Precaución
DA - Forsigtig
SV - Varning
PTBR - Cuidado
TR - Dikkat
EL - Προσοχή
FI - Varoituis
NO - Forsiktig

EN - Importer
DE - Importeur in die EU
NL - Importeur
FR - Importateur
IT - Importatore
ES - Importador
DA - Importør
SV - Importör
PTBR - Importador
TR - İthalatçı
EL - Εισαγωγέας
FI - Maahantuoja
NO - Importør

Rx only

EN - CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician.

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